REAADS[®]
IgA Anti-Phosphatidylserine Semi-Quantitative Test Kit

For *In Vitro* Diagnostic Use

An enzyme-linked immunosorbent assay (ELISA) for the semi-quantitative determination of IgA anti-phosphatidylserine (aPS) antibodies in human serum and citrated plasma (3.2% sodium citrate).

INTENDED USE

Detection and semi-quantitation of IgA anti-phosphatidylserine (aPS) antibodies as an aid for assessing the risk of thrombosis in individuals with systemic lupus erythematosus (SLE) and lupus-like disorders (anti-phospholipid syndrome).

SUMMARY AND EXPLANATION OF THE IGA ANTI-PHOSPHATIDYLSERINE TEST

Anti-phospholipid antibodies are a heterogeneous group of immunoglobulins that bind to several anionic phospholipids, including cardiolipin and phosphatidylserine.^{1,2} High serum levels of anti-phospholipid antibodies are frequently detected in patients with autoimmune (e.g., SLE) and non-autoimmune diseases, as well as in apparently healthy individuals.^{3,4} These antibodies have been associated with an increased risk for recurrent arterial and venous thrombotic events, thrombocytopenia, and fetal loss. These manifestations are the main features of the anti-phospholipid syndrome.^{5,6,7}

Anti-phospholipid antibodies are detected by either ELISAs using cardiolipin as the antigen (anti-cardiolipin antibodies)⁸ or coagulation assays (lupus anticoagulants).⁹ Unlike cardiolipin, phosphatidylserine is a more physiologically relevant phospholipid due to its presence in cell membranes of endothelial cells and platelets¹⁰ and its role in the coagulation cascade.¹¹ The detection of anti-phosphatidylserine (aPS) antibodies by ELISA has been recommended for the serological diagnosis of anti-phospholipid syndrome.^{9,12,13} Patients with positive reactions to both cardiolipin and phosphatidylserine were more likely to have clinical complications than those positive for only one.¹⁴ Higher prevalence and mean serum levels of anti-phosphatidylserine antibodies have been reported in autoimmune patients.¹⁵ In addition, antiphosphatidylserine antibodies in SLE patients correlated with clinical manifestations of anti-phospholipid syndrome^{16,17} and their pathogenic role has been demonstrated in a murine model.¹⁸

Autoimmune anti-phospholipid antibodies require a serum cofactor (β_2 -glycoprotein I) for optimal binding.^{19,20} Antibodies to phosphatidylserine also require β_2 -glycoprotein I.²¹ Bovine serum is provided as the source of cofactor in the formulation of this assay for increased detection of clinically relevant antibodies. The REAADS IgA Anti-phosphatidylserine Test Kit uses the well known ELISA format to detect anti-phosphatidylserine antibodies in human serum or plasma and provides rapid, highly reproducible, accurate, and objective results.

PRINCIPLE OF THE TEST

The test is performed as an indirect ELISA. Diluted serum or citrated plasma samples, calibrator sera, and controls are incubated in phosphatidylserine coated microwells. β_2 -glycoprotein I is provided in the sample diluent. Incubation allows the anti-phosphatidylserine (aPS) antibodies present in the samples to react with the immobilized antigen. After the removal of unbound serum or plasma proteins by washing, antibodies specific for human IgA, labeled with horseradish peroxidase (HRP), are added forming complexes with the phosphatidylserine bound antibodies. Following another washing step, the bound enzyme-antibody conjugate is assayed by the addition of a single solution containing tetramethylbenzidine (TMB) and hydrogen peroxide (H_2O_2) as the chromogenic substrate. Color develops in the wells at an intensity proportional to the serum concentration of IgA anti-phosphatidylserine (aPS) antibodies.

Results are obtained by reading the O.D. (optical density or absorbance) of each well in a spectrophotometer. IgA calibrator sera are provided, expressed as APS (IgA anti-phosphatidylserine) units. The user has the option of running either a single point calibrator or a four-point calibration curve. For single point calibration, dividing the concentration value of the calibrator sera by the O.D. value of the calibrator provides a conversion factor. The O.D. values of the other samples are multiplied by the conversion factor to obtain IgA anti-phosphatidylserine antibody concentrations in APS units. For multipoint calibration, perform a linear regression analysis with calibrator values against calibrator O.D.s. Control and patient results are determined from the calibration curve.

REAGENTS

Store at 2 - 8°C. Do Not Freeze.

Each REAADS IgA anti-Phosphatidylserine 96-microwell Test Kit contains the following reagents (**volumes may vary depending on the kit size and configuration**):

- 12x8 Phosphatidylserine (from porcine brain) coated microwells, with frame.
- 60 mL Sample Diluent I (green solution); contains bovine calf serum*
- 0.250 mL IgA Calibrator Sera* (1-high, 2-moderate, 3-low)(human); see vial label for antibody concentration in APS units. Calibrator 3 should be used when performing single point calibration.
- 0.250 mL IgA Positive Control Serum* (human); see vial label for expected APS range.
- 0.250 mL Normal Control Serum* (human); see vial label for expected APS range.
- 15 mL anti-human IgA (goat) HRP-Conjugated Antibody Solution (orange solution); contains 0.01% thimerosal and gentamycin sulfate as preservatives.
- 15 mL One Component Substrate (TMB/ H_2O_2); ready to use.
- 15 mL Stopping Solution (0.36 N sulfuric acid)
- 2 bottles (30 mL) Wash Concentrate (33X PBS)

***CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use

1. Human source material used to prepare the calibrators and controls included in this kit has been tested and shown to be negative for antibodies to HBsAg, HCV, and HIV I & 2 by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Certain components of this product contain sodium azide as a preservative. Sodium azide has been reported to form lead and copper azides when left in contact with these metals. These metal azides are explosive. Any solutions containing azide must be thoroughly flushed with copious amounts of water to prevent the build-up of explosive metal azides in the plumbing system.

6. One component substrate can cause irritation to the eyes and skin. Absorption through the skin is possible. Use gloves when handling substrate and wash thoroughly after handling. Keep reagent away from ignition sources. Avoid contact with oxidizing agents.
7. Certain components are labeled with the following:
Irritating to eyes (R 36). Avoid contact with skin (S 24). Avoid contact with eyes (S 25). In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S 26). If swallowed, seek medical advice immediately and show this container or label (S 46).

Irritant . Biological Risk .

SPECIMEN COLLECTION AND PREPARATION

Serum is the preferred sample matrix. Blood should be collected by venipuncture, and the serum separated from the cells by centrifugation after clot formation. If not tested immediately, specimens should be stored at 2 to 8°C. If specimens are to be stored for more than 72 hours, they should be frozen at -20°C or below. Avoid repeated freezing and thawing. Do not use hemolyzed, icteric, or lipemic serum or plasma as these conditions may cause aberrant results. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

Plasma collected with 3.2% sodium citrate may be used. Blood should be collected by venipuncture and the plasma immediately separated from the cells by centrifugation at 1500g for 10 minutes. The supernatant must be carefully removed after centrifugation to avoid contamination with platelets. Repeating the centrifugation and separation steps may be advisable to minimize platelet contamination. Lysed or aged platelets can react with anti-phospholipid antibodies leading to aberrant results. If not tested immediately, plasma samples should be stored as described for serum.

INSTRUCTIONS FOR USE

Materials Provided:

REAADS IgA anti-Phosphatidylserine Test Kit; see “Reagents” for a complete list.

Materials Required but not Supplied

- Reagent grade water to prepare PBS wash solution and to zero or blank the plate reader during the final assay step
- Graduated cylinders
- Precision pipettors capable of delivering between 5 µL and 1000 µL, with appropriate tips
- Miscellaneous glassware appropriate for small volume handling
- Flask or bottle, 1 liter
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semi-automated washing system
- Disposable gloves
- Plate reading spectrophotometer capable of reading absorbance at 450 nm (with a 650 nm reference if available)
- Multichannel pipettors capable of delivering to 8 wells simultaneously

Procedural Notes

1. Bring serum samples and kit reagents to room temperature (18-26°C) and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage as soon as possible.
2. All dilutions of calibrators, controls, and test sera must be made just prior to use in the assay.
3. A single water blank well can be set up on each plate with each run. No sample or kit reagents are to be added to this well. Instead, add 200 µL of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to zero or blank against an air or a water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution into the bottom of the microwells from a plastic squeeze bottle with a wide tip. Wash solution in the water blank well will not interfere with the procedure. An automated microtiter plate washing system can also be used.
5. **IMPORTANT:** Failure to adequately remove residual PBS can cause inconsistent color development of the Substrate Solution.

6. Use a multichannel pipettor capable of delivering to 8 wells simultaneously when possible. This speeds the process and provides more uniform incubation and reaction times for all wells.
7. Careful controlled timing of all steps is critical. All calibrators, controls, and samples must be added within a five minute period. Batch size of samples should not be larger than the amount that can be added within this time period.
8. For all incubations, the start of the incubation period begins with the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below normal room temperature (18-26°C) may contribute to inaccurate results.
11. Avoid contamination of reagents when opening and removing aliquots from the primary vials.
12. Do not use Tween 20 or other detergents in this assay.
13. Do not use kit components beyond expiration date.
14. Do not use kit components from different kit lot numbers.

REAGENT PREPARATION

Wash Solution (PBS): Measure 30 mL of Wash Concentrate (33X PBS) and dilute to 1 L with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused PBS solution in the refrigerator at 2-8°C. Discard if the solution shows signs of microbial or cross-contamination.

ASSAY PROCEDURE

1. The assay can be performed with a single point calibration (Calibrator 3) or a four-point calibration curve (Calibrators 1, 2, and 3 plus sample diluent/reagent blank as Calibrator 4 equal to 0 APS units). One well should be run as a reagent blank control with the single point and multipoint calibration method; Sample Diluent I without serum or plasma is added to the well. This well will be treated the same as sample wells in subsequent assay steps.
2. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
3. Prepare a 1:50 dilution of the calibrators, controls, and patient samples in sample diluent (green solution); e.g., 10 μ L sample added to 490 μ L Sample Diluent I equals a 1:50 sample dilution.
4. Add 100 μ L of diluted calibrators (including the reagent blank/Calibrator 4), controls, and patient samples to the appropriate microwells.
5. Incubate 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and empty the sample fluid. Do not allow samples to contaminate other microwells.
6. Wash 4 times with PBS. Each well should be filled with PBS per wash. Invert microwells between each wash to empty fluid. Use a snapping motion of the wrist to shake the liquid from the wells. The frame must be squeezed at the center on the top and bottom to retain microwell modules during washing. Blot on absorbent paper to remove residual wash fluid. Do not allow wells to dry out between steps.
7. Add 100 μ L IgA Conjugate to the calibrators, controls, reagent blank, and patient samples wells.
8. Incubate for 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and dump the conjugate solution.
9. Wash 4 times with PBS as in step 7. Use a snapping motion to drain the liquid and blot on absorbent paper after the final wash. Do not allow the wells to dry out.
10. Add 100 μ L Substrate to each well and incubate for 10 minutes at room temperature. Add the substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
11. Add 100 μ L Stopping Solution (0.36 N sulfuric acid) to each well to stop the enzyme reaction. Be sure to add the Stopping Solution to the wells in the same order and at the same rate as the Substrate was added. Blue Substrate will turn yellow and colorless solution will remain colorless. Blank or zero the plate reader against an air or a water blank well. Read the O.D. of each well at 450 nm (and 650 nm reference if dual beam). The O.D. values should be measured within 5 minutes after the addition of Stopping Solution.

RESULTS

Single Point Calibration


1. Calculate the mean O.D. values if duplicates of Calibrator 3, controls, and patient samples were performed.
2. Divide the concentration value of Calibrator 3 (printed on the vial label) by the O.D. or mean O.D. value of the calibrator serum to obtain the conversion factor.
3. Multiply the O.D. or mean O.D. values for each of the controls and patient samples by the conversion factor to obtain a concentration value in APS units.

$$\text{Conversion Factor} = \frac{\text{IgA Anti-Phosphatidylserine Concentration of Calibrator 3 (APS units)}}{\text{Absorbance Value of the Calibrator 3 (O.D.)}}$$
$$\text{Anti-phosphatidylserine concentration of sample} = \text{Conversion Factor} \times \text{Absorbance of the sample (O.D.)}$$

4. The Conversion Factor must be calculated for each assay run. Using a Conversion Factor from another assay will invalidate the results.

Multi-Point Curve Calibration

1. Calculate the mean O.D. values if duplicates of the calibrators, controls and patient samples were performed.
2. Perform linear regression analysis with the four calibrator values (See vial labels for APS units. Calibrator 4 [Sample Diluent] is equal to 0 APS units) against the O.D.s for each calibrator.
3. The calibrator curve can be plotted either automatically using a validated software program or manually with graph paper. It is recommended to use a zero intercept when generating the regression line to avoid negative values. If this option is not available, any negative values should be reported as zero units. When generating the curve manually, draw a best fit line through the plotted points using a zero intercept.
4. Determine the control and patient sample values from the calibrator curve.
5. Example of a multi-point curve calibration.

Using the example calibration curve provided, a specimen O.D. of 1.000 at 450 nm would correspond to a calculated value of 26.2 units. The calibration curve provided is an example only and should not be used to calculate patient results. A new calibration curve should be performed with every test run.

QUALITY CONTROL

1. The O.D. value of the Calibrator 3 should be at least 0.400 to assure that the kit is functioning properly. Calibrator 3 O.D. readings of less than 0.400 may indicate that the kit is no longer suitable for use.
2. The O.D. of the Calibrator 4 or reagent blank should be less than 0.150 when the spectrophotometer has been blanked against the water well. Readings greater than 0.150 may indicate possible reagent contamination or inadequate plate washing.
3. The anti-phosphatidylserine antibody values obtained for the control sera should be within the ranges indicated on the container labels. Occasional small deviations outside these ranges are acceptable.
4. O.D. values for duplicates of the controls or patient samples should be within 20% of the mean O.D. value for samples with absorbance readings greater than 0.200.
5. Each laboratory should periodically determine their own normal cut-off values for the appropriate population of patients. See Performance Characteristics, Clinical Specificity, as an example.
6. Samples with IgA anti-phosphatidylserine values greater than 60 APS units may be reported as "greater than 60 APS".
7. Assure that all quality control parameters have been met (see Quality Control section) before reporting test results.

NORMAL RANGE

Serum samples from 100 healthy blood donors were tested for IgA anti-phosphatidylserine (aPS) antibodies. The normal range was established as less than 20 APS units (mean + 2 SD).

PERFORMANCE CHARACTERISTICS

Normal Samples:

Serum samples from 100 healthy blood donors were assayed in three lots for the presence of IgA anti-phosphatidylserine antibodies. Using the established cutoff value of 20 APS units, this normal population demonstrated a mean 95% specificity (mean value 11.5 APS units).

Serum samples from 42 infectious disease (Syphilis) patients, 42 Progressive Systemic Sclerosis (PSS) patients, and 42 Rheumatoid Arthritis (RA) patients were assayed for the presence of IgA anti-phosphatidylserine antibodies. These groups demonstrated mean values of 13.5, 15.8, and 14.7 APS units respectively. Results of these groups are summarized in the table below.

Summary	Healthy Blood Donors	Syphilis	PSS	RA
# of Samples (n)	100	42	42	42
Mean Value (APS units) + SD	11.5 + 5.0	13.5 + 8.6	15.8 + 7.9	14.7 + 63
% Negative (Specificity)	95%	88.1%*	85.7%**	81.0%**

*p>0.05 versus healthy blood donors

**p<0.05 versus healthy blood donors

Clinical Sensitivity:

Systemic Lupus Erythematosus (SLE)

Serum samples from 34 unselected (consecutive) patients with SLE were tested with the kit. Of the samples four (sensitivity 11.7%) were positive for IgA aPS (mean value 12.6 APS units). A good correlation ($r=0.976$) was seen between IgA anti-Phosphatidylserine and IgA anti-Cardiolipin levels in this group.

Serum samples from 15 selected female patients with SLE who had a clinical history of thrombosis, thrombocytopenia, or recurrent fetal loss (Secondary APS) were evaluated for antibodies to IgA anti-Phosphatidylserine. Six of these samples (sensitivity 40%) were positive in this population (mean value 19.5 APS units). A good correlation ($r=0.985$) was seen between IgA anti-Phosphatidylserine and IgA anti-Cardiolipin levels in this group.

SLE Controls

Serum samples from nine selected female patients with SLE who were not known to have had a thrombotic event, nor any other feature of the anti-phospholipid syndrome (controls), were evaluated for antibodies to IgA anti-Phosphatidylserine. None of the samples were found positive for IgA anti-phosphatidylserine antibodies, with a mean value of 11.2 APS units. This group showed good correlation ($r=0.986$) to IgA anti-cardiolipin antibody testing.

Primary Anti-Phospholipid Syndrome

Serum samples from 11 patients with the diagnosis of primary antiphospholipid syndrome (APS) were tested on the IgA anti-Phosphatidylserine assay. Four of these samples (36.4%) were found positive for anti-phosphatidylserine antibodies with a mean value of 21.1 APS units. This group showed good correlation ($r=0.932$) to IgA anti-Cardiolipin antibody testing. A summary of the sensitivity testing is presented below.

Summary	Selected SLE			
	Unselected SLE	SLE Controls	SecondaryAPS (SLE)	PrimaryAPS
# of Samples (n)	34	9	15	11
Mean Value (APS units) + SD	12.6 + 7.3	11.2 + 4.5	19.5 + 10.8	21.1 + 15
% Positive (Sensitivity)	11.7%	0%*	40%**	36.4%**

* $p>0.05$ versus healthy blood donors

** $p<0.05$ versus healthy blood donors

Precision:

Two samples with known IgA aPS values (one low and one high) were assayed in 23 replicates on three different occasions. The mean intraassay and interassay coefficients of variation (CVs) are presented in the following table. The reported intraassay coefficient of variation is the mean of the three separate intraassay CVs. Interassay CV is the coefficient of variation obtained from three plates from one lot.

Sample IgA aPS Value	Intraassay CV	Mean Interassay CV
Low (< 10 APS units)	11.3%	7.0%
High (> 30 APS units)	4.7%	4.6%

LIMITATIONS OF THE TEST

The IgA anti-phosphatidylserine antibody concentration values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures. If clinical findings suggest the presence of anti-phospholipid antibodies and the patient is negative for antiphosphatidylserine antibodies, some investigators recommend testing for anti-cardiolipin and anti- β_2 GPI antibodies and the lupus anticoagulant to confirm the negative result. A patient is considered positive for anti-phospholipid antibodies if one or all of the tests give positive results .

Anti-phosphatidylserine antibodies can appear transiently at low levels during many infections.^{22,23} If a patient first tests positive while there are clinical signs of infection, the test should be repeated after an interval of six months.

WARRANTY

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone (303)-457-4345, fax (303)-457-4519, email: techsupport@corgenix.com, or contact a Corgenix authorized distributor.

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REFERENCES

1. McNiel HP, Chesterman CN, Krilis SA. Immunology and clinical importance of antiphospholipid antibodies. *Adv Immunol* 1991; 49:193-280.
2. Harris EN. Annotation. Antiphospholipid antibodies. *Br J Haematol* 1990; 74:1-9.
3. Love PE, Santoro SA. Antiphospholipid antibodies: Anticardiolipin and the lupus anticoagulant in systemic lupus erythematosus (SLE) and in non-SLE disorders. Prevalence and clinical significance. *Ann Intern Med* 1990; 112:682-698.
4. Sammaritano LR, Gharavi AE, Lockshin MD. Antiphospholipid antibody syndrome: immunologic and clinical aspects. *Seminars in Arthritis and Rheumatism* 1990; 20:81-96.
5. Hughes GRV, Harris EN, Gharavi AE. The anticardiolipin syndrome (Editorial). *J Rheumatol* 1986; 13:486-489.
6. Hughes GRV. The antiphospholipid syndrome: ten years on. *Lancet* 1993; 342:341-344.
7. Wilson, WA, Gharavi AE, Koike T, et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. *Arthritis Rheum* 1999;42:1309-1311.
8. Loizou S, McCrea JD, Rudge AC, et al. Measurement of anti-cardiolipin antibodies by an enzyme linked immunosorbent assay (ELISA): standardization and quantitation of results. *Clin Exp Immunol* 1985; 62:738-745.
9. Triplett DA. Assays for detection of antiphospholipid antibodies. *Lupus* 1994; 3:281-287.
10. Schick PK, Kurica KB, Chacko GK. Location of phosphatidylethanolamine and phosphatidylserine in the human platelet plasma membrane. *J Clin Invest* 1976; 57:1221-1226.
11. Branch DW, Rote NS, Dostal DA, Scott JR. Association of lupus anticoagulant with antibody against phosphatidylserine. *Clin Immunol Immunopathol* 1987; 42:63-75.
12. Wilson WA, Faghiri Z, Taheri F, Gharavi AE. Significance of IgA antiphospholipid antibodies. *Lupus* 1998; 7:Suppl 2, S110-S113.
13. Bick RL, Baker WF. Antiphospholipid syndrome and thrombosis. *Seminars in Thrombosis and Hemostasis* 1999; 25: 333-350.
14. Barna LK, Triplett DA, Foster WM, Gaddis ML. A study of relationships among anti-phosphatidylserine and anticardiolipin antibodies, the lupus anticoagulant and clinical complications. *Lupus* 1994; 3:357.
15. Keedy K, Santos M, Lopez L. Prevalence of antibodies to anti-phosphatidylserine in autoimmune diseases. *Lupus* 1994; 3:344.
16. Keedy K, Santos M, Lopez L. Prevalence and clinical significance of anti-phosphatidylserine antibodies. *Arthritis Rheum* 1994; 37:S392.
17. Jansen NL, Snell JR, Moy JN. Myocardial infarction as the presenting manifestation of systemic lupus erythematosus with antiphosphatidylserine antibodies. *Annals of Allergy, Asthma, and Immunology* 1996; 76:266-268.
18. Blank M, Tincani A, Shoenfeld Y. Induction of experimental antiphospholipid syndrome in naive mice with purified IgG anti-phosphatidylserine antibodies. *J Rheumatol* 1994; 21:100-104.
19. McNiel H, Simpson R, Chesterman C, Krilis S. Antiphospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β 2-glycoprotein I (apolipoprotein H). *Proc Natl Acad Sci USA* 1990; 87:4120.
20. Kandiah D, Krilis S. Beta2-glycoprotein I. *Lupus* 1994; 3:207-212.
21. Keedy K, Santos M, Lopez L. Anti-phosphatidylserine antibodies require β 2-glycoprotein I as cofactor in ELISA. *Lupus* 1994; 3:327.
22. Vaarala O, Palosuo T, Kleemola M, Aho K. Anticardiolipin response in acute infections. *Clinical Immunology and Immunopathology* 1986; 41: 8-15.
23. Maclean C, Flegg PJ, Kilpatrick DC. Anti-cardiolipin antibodies and HIV infection. *Clinical Experimental Immunology*. 1990; 81: 263-266.

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Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020, USA

REAADS® is a registered trademark of Corgenix, Inc.
 © 2004, Corgenix, Inc.

10207 08
 Effective:2006-10-12

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MT Promedt Consulting GmbH