LABELING

Title: **REAADS IgM ANTI-PROTHROMBIN
PACKAGE INSERT**

Document Number: 10294

Supersedes/Date: 2007-04-12

CO No.: 07-0630

Revision: 08

QA Approval/Date: K.Hassler 2007-11-07

Effective Date: 2008-01-11

Dimensions: Flat 8½" wide X 11" long

Folded / Booklet Format 5 ½" X 8 ½"

Colors: Black Text

Paper Stock: White

REAADS[®]
IgM Anti-Prothrombin Semi-Quantitative Test Kit

For *In Vitro* Diagnostic Use

An enzyme-linked immunosorbent assay (ELISA) for the semiquantitative determination of IgM anti-prothrombin (aPT) antibodies in human serum or citrated plasma (3.2% sodium citrate).

INTENDED USE

For the detection and semi-quantitation of IgM antiprothrombin (aPT) antibodies in individuals with systemic lupus erythematosus (SLE) and lupus-like disorders (e.g., antiphospholipid syndrome).

SUMMARY AND EXPLANATION OF THE ANTI-PROTHROMBIN TEST

Antiphospholipid antibodies are a heterogeneous group of immunoglobulins (IgG, IgM, IgA) that bind to several anionic phospholipids (e.g., cardiolipin, phosphatidylserine), to phospholipid-protein complexes, and to certain proteins (cofactors) in the absence of anionic phospholipids.¹⁻³ Several plasma proteins, mostly associated with the coagulation system and with strong phospholipid binding properties, have been identified as antiphospholipid cofactors. Beta 2 glycoprotein I (β_2 GPI) and prothrombin are the most common and extensively studied antiphospholipid cofactors.⁴⁻⁶ Most autoimmune antiphospholipid antibodies, which are strongly associated with thrombosis, require a cofactor for optimal interaction with the antigen. The epitope responsible for this interaction is localized on the cofactor molecule.^{7,8} Anti- β_2 GPI antibodies have been shown to be more specific for thrombosis than anti-cardiolipin (aCL) antibodies.⁹⁻¹¹ Many clinical laboratories now include the detection of cofactor antibodies (e.g., anti- β_2 GPI) as part of their diagnostic panel for the evaluation of antiphospholipid antibodies. More recently, antibodies to the other protein cofactor (prothrombin) have been reported in patients with the antiphospholipid syndrome (APS)¹² and have been proposed to be included in the serologic evaluation of antiphospholipid antibodies. Clinically, elevated serum levels of these antibodies are associated with an increased risk for APS, characterized by recurrent arterial or venous thrombosis, thrombocytopenia, and/or fetal loss. APS may be found in patients with autoimmune diseases (e.g., SLE or secondary APS) as well as in individuals without an obvious underlying disease (primary APS).^{1,13,14}

Antiphospholipid antibodies are classified according to the laboratory method used for their detection. Lupus anticoagulants (LA) are a group of antiphospholipid antibodies that inhibit phospholipid-dependent coagulation assays. LAs have been shown to correlate better with thrombosis (APS) than classic antiphospholipid ELISA assays (e.g., anticardiolipin and anti-phosphatidylserine).¹⁵ LAs also require protein cofactors (β_2 GPI or prothrombin) for optimal immunologic binding.¹⁶ The majority of LA samples contain both antibodies (anti- β_2 GPI and anti-prothrombin) which have been identified as responsible for the LA activity.¹⁷ However, it is not infrequent to find LA positive samples that contain only anti-prothrombin (aPT) antibodies. Galli and collaborators^{18,19} described different coagulation profiles for samples containing anti- β_2 GPI or aPT antibodies. LA activity defined by an abnormal kaolin clotting time (KCT) seems to be caused by aPT antibodies while LA activity defined by an abnormal dRVVT clotting time seems to be caused by anti- β_2 GPI antibodies. Furthermore, patients with these abnormal clotting profiles showed different incidence of thrombosis.²⁰ The exact pathogenic role of aPT antibodies in thromboembolic disease and as a serologic marker for assessing APS is now under intensive investigation.²¹⁻²³

The REAADS anti-prothrombin (aPT) ELISA test kit uses purified human prothrombin as antigen in a well-known ELISA format to detect IgM anti-prothrombin antibodies in human serum or citrated plasma in the absence of other exogenous cofactors or phospholipids. High serum or plasma levels of aPT antibodies may add valuable information in the laboratory assessment of antiphospholipid antibodies.

PRINCIPLE OF THE TEST

The test is performed as an indirect ELISA. Diluted serum or citrated plasma samples, calibrator, and controls are incubated in microwells coated with purified human prothrombin. Incubation allows the anti-prothrombin antibodies present in the samples to react with the immobilized antigen. After the removal of unbound proteins by washing, antibodies specific for human IgM labeled with horseradish peroxidase (HRP) are added forming complexes with the prothrombin bound antibodies. Following another washing step, the bound enzyme-antibody conjugate is assayed by the addition of a single solution containing tetramethylbenzidine (TMB) and hydrogen peroxide (H₂O₂) as the chromogenic substrate. Color develops in the wells at an intensity proportional to the concentration of anti-prothrombin antibodies.

Results are obtained by reading the optical density (O.D. or absorbance) of each well in a spectrophotometer. Calibrator sera are provided with the IgM anti-prothrombin antibody concentration expressed in M units. The user has the option of running either a single point calibrator or a four-point calibration curve. For single point calibration, dividing the concentration value of the calibrator sera by the O.D. value of the calibrator provides a conversion factor. The O.D. values of the other samples are multiplied by the conversion factor to obtain IgM anti-prothrombin antibody concentrations in M units. For multipoint calibration, perform a linear regression analysis with calibrator values against calibrator O.D.'s. Controls and patient results are determined from the calibration curve.

REAGENTS

Store at 2 - 8°C. Do Not Freeze.

Each REAADS IgM anti-Prothrombin 96-microwell Test Kit contains the following reagents (volumes may vary depending on the kit size and configuration):

- 12 x 8 Stabilized Prothrombin (human) coated microwells with Frame
- 60 mL Sample Diluent II (blue-green solution)
- 3 vials (0.25 mL each) IgM anti-Prothrombin Calibrator (1- high, 2-moderate, 3-low)(human); see vial label for antibody concentration in M units.* Calibrator 2 should be used when performing single point calibration.
- 0.25 mL IgM anti-Prothrombin Positive Control (human); see vial label for expected M unit range*
- 0.25 mL Normal Control (human); see vial label for expected M unit range*
- 15 mL anti-Human IgM (goat) HRP-Conjugated Antibody Solution (red solution)
- 15 mL One-Component Substrate Solution; contains 3,3',5,5'-tetramethylbenzidine, 0.01% hydrogen peroxide in acidic buffer
- 15 mL Stopping Solution (0.36 N sulfuric acid)
- 2 bottles (30 mL) Wash Concentrate (33X PBS/Tween)

***CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use

1. Human source materials used to prepare the calibrator and controls included in this kit have been tested and shown to be negative for antibodies to HBsAg, HCV and HIV 1 & 2 by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Certain components of this product contain sodium azide as a preservative. Sodium azide has been reported to form lead and copper azides when left in contact with these metals. These metal azides are explosive. Any solutions containing azide must be thoroughly flushed with copious amounts of water to prevent the build-up of explosive metal azides in the plumbing system.

6. One-component substrate can cause irritation to the eyes and skin. Use gloves, eye protection, and lab coat when handling. Wash hands thoroughly after handling. Keep reagent away from ignition sources. Avoid contact with oxidizing agents.
7. Certain components are labeled with the following:
Irritating to eyes (R 36). Avoid contact with skin (S 24). Avoid contact with eyes (S 25). In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S 26). If swallowed, seek medical advice immediately and show this container or label (S 46).

Irritant  Biological Risk .

SPECIMEN COLLECTION AND PREPARATION

Serum or citrated plasma (3.2% sodium citrate) should be used as the sample matrix. Blood should be collected by venipuncture and the serum separated from the cells by centrifugation after clot formation.

If citrated plasma is to be used, blood should be collected by venipuncture and the plasma separated from the cells immediately by centrifugation at 1500g for 10 minutes. The supernatant must be carefully removed after centrifugation to avoid contamination with platelets. Repeating the centrifugation and separation steps may be advisable to minimize platelet contamination. Lysed or aged platelets can lead to aberrant results.

If not tested immediately, the specimen should be stored at 2 to 8°C. If specimens are to be stored for more than 72 hours, they should be frozen at -20°C or below. Avoid repeated freezing and thawing. Do not use hemolyzed, icteric, or lipemic serum or plasma as these conditions may cause aberrant results. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

INSTRUCTIONS FOR USE

Materials Provided:

REAADS IgM anti-Prothrombin Test Kit; see “Reagents,” section for a complete list.

Materials Required but not Supplied:

- Reagent grade water to prepare wash solution (2L) and to zero or blank the plate reader during the final assay step
- Graduated cylinders
- Precision pipettors capable of delivering between 5 µL and 1000 µL, with appropriate tips
- Miscellaneous glassware appropriate for handling small volumes
- Flasks or bottles, 1 L
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semiautomated washing system
- Disposable gloves
- Plate-reading spectrophotometer capable of reading absorbance at 450 nm (650 nm reference if dual beam)
- Multichannel pipettors capable of delivering to 8 wells simultaneously
- Microdilution tubes and a 96-well microdilution tube holder for sample dilutions and rapid delivery to microwell plate

Procedural Notes

1. Bring serum or plasma samples and kit reagents to room temperature and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage as soon as possible.
2. All dilutions of the calibrator, controls, and test sera or plasma must be made just prior to use in the assay.
3. A single water blank well can be included in each plate with each run. No sample or kit reagents are to be added to this well. Instead, add 200 µL of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to “zero” or “blank” against air or a water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution from a plastic squeeze bottle with a wide tip into the bottom of the microwells. Wash solution in the water blank well will not interfere with the procedure. An automated plate washing system can also be used.

5. **IMPORTANT:** Failure to adequately remove residual wash solution can cause inconsistent color development of the Substrate Solution.
6. Use a multichannel pipettor capable of delivering 100 μ L each to 8 wells simultaneously when possible. This speeds the process and provides more uniform incubation and reaction times for all wells.
7. Carefully controlled timing of all steps is critical. All calibrators, controls, and samples must be added within a five-minute period. Batch size of samples should not be larger than the amount that can be added within this time period.
8. For all incubations, the start of the incubation period begins with the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below normal room temperature (18° to 26°C) may contribute to inaccurate results.
11. Avoid microbial and cross-contamination of reagents when opening and removing aliquots from the primary vials.
12. Do not use kit components beyond expiration date.
13. Do not mix or use kit components from different kit lot numbers.

REAGENT PREPARATION

Wash Solution (PBS/Tween): Measure 30 mL of wash concentrate and dilute to 1 L with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused wash solution in the refrigerator at 2-8°C. Discard if the solution shows signs of microbial or cross-contamination.

Assay Procedure

1. The assay can be performed with a single point calibration (Calibrator 2) or a four-point calibration curve (Calibrators 1, 2, and 3 plus sample diluent/reagent blank as Calibrator 4 equal to 0 M units). A reagent blank control should also be run with the single point calibration method. Sample Diluent without serum or plasma is added to the well. This well will be treated the same as a control or patient sample in subsequent assay steps.
2. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
3. Prepare a 1:51 dilution of the calibrators, controls, and patient samples in sample diluent (blue-green solution); e.g., 10 μ L sample added to 500 μ L sample diluent equals a 1:51 sample dilution. Mix well.
4. Add 100 μ L of diluted calibrators (including the reagent blank/Calibrator 4), controls, and patient sample(s) to the appropriate microwells.
5. Incubate 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and dump the sample fluid. Do not allow samples to contaminate other microwells. To retain microwells during inversion and washing, squeeze the frame along the longer sides.
6. Wash 4 times with wash solution. Completely fill each well with wash solution. Invert the microwells to empty the fluid. Use a snapping motion of the wrist to shake the liquid from the wells. Repeat the process three additional times. Invert the wells on absorbent paper and tap to blot residual wash solution. Do not allow the wells to dry out at any time.
7. Add 100 μ L anti-Human IgM HRP-Conjugated Antibody Solution (red) to the wells.
8. Incubate for 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and dump the conjugate solution.
9. Wash 4 times with wash solution, as described in step 6. Do not allow the wells to dry out.
10. Add 100 μ L One-Component Substrate Solution to each well and incubate for 10 minutes at room temperature. Add substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
11. Add 100 mL Stopping Solution (0.36 N sulfuric acid) to each well to stop the enzyme reaction. Be sure to add the acid to the wells in the same order and at the same rate as the Substrate was added. Blue Substrate will turn yellow and colorless solution will remain colorless. Blank or zero the plate reader against air or a water blank well. Read the O.D. of each well at 450 nm (and 650 nm reference if dual beam). The O.D. values should be measured within 5 minutes of the addition of the Stopping Solution.

Results

Single Point Calibration

1. Calculate the mean O.D. reading if duplicates of Calibrator 2, Controls, and patient samples were performed.
2. Divide the concentration value of Calibrator 2 (printed on vial label) by its O.D. or mean O.D. reading to obtain the conversion factor.
3. Multiply the O.D. or mean O.D. reading for each of the controls and patient samples by the conversion factor to obtain an anti-prothrombin (aPT) antibody concentration value expressed in M units.


Conversion Factor	=	$\frac{\text{IgM aPT Concentration of Calibrator 2 (in M units)}}{\text{O.D. or Mean O.D. Reading of the Calibrator 2}}$
IgM aPT Concentration of Sample	=	Conversion Factor X O.D. or Mean O.D. of the Sample

4. The Conversion Factor must be calculated for each assay run. Using a Conversion Factor from another assay will invalidate the results.

MULTI-POINT CURVE CALIBRATION

1. Calculate the mean O.D. values if duplicates of the calibrators, controls and patient samples were performed.
2. Perform linear regression analysis with the four calibrator values (See vial labels for M units. Calibrator 4 [sample diluent] is equal to 0 M units) against the O.D.'s for each calibrator.
3. The calibrator curve can be plotted either automatically using a validated software program or manually with graph paper. It is recommended to use a zero intercept when generating the regression line to avoid negative values. If this option is not available, any negative values should be reported as zero units. When generating the curve manually, draw a best fit line through the plotted points using a zero intercept.
4. Determine the control and patient sample values from the calibrator curve.

For Example:

Using the example calibration curve provided, a specimen O.D. of 0.860 at 450 nm would correspond to a calculated value of 26.2 units. The calibration curve provided is an example only and should not be used to calculate patient results. A new calibration curve should be performed with every test run.

Quality Control

1. The O.D. reading of the Calibrator 2 should be ≥ 0.600 to assure that the kit is functioning properly. Calibrator 2 O.D. readings of less than 0.600 may indicate that the kit is no longer suitable for use.
2. The O.D. of Calibrator 4 or reagent blank should be less than 0.100 when the spectrophotometer has been blanked against the water well. Readings greater than 0.100 may indicate possible reagent contamination or inadequate plate washing.
3. The anti-prothrombin (aPT) values obtained for the control samples should be within the ranges indicated on the vial labels. Occasional small deviations outside these ranges are acceptable.
4. O.D. readings for duplicate wells of the controls or patient samples should be within 20% of the mean O.D. reading for that sample if the absorbance readings are greater than 0.200.
5. Each laboratory should determine its own normal cut-off values for the appropriate population of patients.
6. Samples with anti-prothrombin (aPT) values greater than 100 M units may be reported as “greater than 100 M units.”
7. Ensure that all quality control parameters have been met before reporting test results.

NORMAL RANGE

Serum and citrated plasma (3.2% sodium citrate) samples from 100 healthy blood donors were tested for IgM antiprothrombin (aPT) antibodies. The following normal range was established:

- Less than 20 M units

PERFORMANCE CHARACTERISTICS

Clinical Specificity

Normal Serum Samples:

Serum samples from 100 healthy blood donors were assayed for the presence of IgM anti-prothrombin (aPT) antibodies on three lots of REAADS IgM aPT ELISA test kit. Using the established cut-off value of 20 M units, this normal serum population demonstrated 97% specificity (mean value = 7.9 M units, SD = 6.0).

Normal Plasma Samples:

Citrated plasma (3.2% sodium citrate) samples from 100 healthy blood donors were assayed for the presence of IgM anti-prothrombin (aPT) antibodies on two lots of REAADS IgM aPT ELISA test kit. Using the established cut-off value of 20 M units, this normal plasma population demonstrated 97% specificity (mean value = 6.3 M units, SD = 4.3).

Diseased Serum Samples:

Serum samples from 42 infectious disease (syphilis), 10 Osteoarthritis (OA), 42 progressive systemic sclerosis (PSS), and 42 rheumatoid arthritis (RA) patients were assayed for the presence of IgM anti-prothrombin (aPT) antibodies on REAADS IgM aPT ELISA test kit. These patient groups demonstrated statistically similar results compared to the healthy blood donor population (mean values = 10.3, 8.4, 11.2, and 10 M units respectively). Results are summarized in the table below.

Serum Samples	Healthy	Infectious (Syphilis)	Osteo-arthritis (OA)	PSS	RA
# of Samples (n)	100	42	10	42	42
Mean Value \pm SD (M units)	7.9 \pm 6.0	10.3 \pm 8.3	8.4 \pm 5.4	11.2 \pm 7.8	10 \pm 7.1
% Negative (Specificity)	97%	88.1%	90%	90.5%	95.2%

Clinical Sensitivity

Autoimmune Disease:

Serum samples from a group of 40 patients with autoimmune diseases were selected for their positive reactivity to antiphospholipid antibodies and included some patients with SLE and clinical manifestations suggestive of the antiphospholipid syndrome. These serum samples were tested for the presence of IgM aPT antibodies on REAADS IgM aPT ELISA test kit. Using the cut-off of 20 M units, 45% of the samples reacted positive for IgM aPT antibodies (mean value = 23.5 M units, SD = 19.3). The mean IgM aPT value of this population was statistically different from the healthy blood donors ($p < 0.001$).

Systemic Lupus Erythematosus (SLE):

Serum samples from 41 unselected (consecutive) patients with SLE were tested for the presence of IgM aPT antibodies on the REAADS IgM aPT ELISA test kit. Using the cut-off of 20 M units, 12.2% of the samples reacted positive for IgM aPT antibodies (mean value = 11.8 M units, SD = 9.4). The mean IgM aPT value of this population was statistically different from the healthy blood donors ($p = 0.0044$).

Primary Antiphospholipid Syndrome (APS):

Serum samples from 11 patients with primary antiphospholipid syndrome (APS) were tested for the presence of IgM aPT antibodies on REAADS IgM aPT ELISA test kit. Using the cutoff of 20 M units, 27.3% of the samples reacted positive for IgM aPT antibodies (mean values = 15.7 M units, SD = 12.6). The mean IgM aPT value of this population was statistically different from the healthy blood donors ($p < 0.001$).

Serum Samples	Autoimmune	Unselected SLE	Primary APS
# of Samples (n)	40	41	11
Mean Value \pm SD (M units)	23.5 \pm 19.3	11.8 \pm 9.4	15.7 \pm 12.6
% Positive (Sensitivity)	45%	12.2%	27.3%

Lupus Anticoagulant (LA):

Sixty-one (61) citrated plasma (3.2% sodium citrate) samples with LA activity from three clinical institutions were tested for the presence of IgM aPT antibodies on REAADS IgM aPT ELISA test kit. The determination of LA activity was performed by each institution following published criteria for the diagnosis of lupus anticoagulant.²⁴ Fifty-six (56) plasma samples with abnormal coagulation profiles but not meeting the diagnostic criteria for LA (controls) were also tested for the presence of IgM aPT antibodies. Using the cut-off of 20 M units, 24.6% of the samples of the LA positive tested positive for IgM aPT antibodies while only 5.4% of the controls tested positive. The mean IgM aPT value of the LA positive group was statistically different from the controls ($p = 0.001$). A summary of LA sensitivity testing is presented below.

Plasma Samples	Healthy	LA Controls	LA Positive
# of Samples (n)	100	56	61
Mean Value \pm SD (M units)	6.3 \pm 4.3	6.9 \pm 6.7	17.8 \pm 23.9
% Positive (Sensitivity)	3%	5.4%	24.6%

Comparison with antiphospholipid status (APL) of the sample

Healthy and disease state populations were tested for IgM isotypes of aPT, anticardiolipin (aCL), antiphosphatidylserine (aPS) and anti-Beta-2 Glycoprotein I ($\alpha\beta_2$ GPI) antibodies²⁵. Samples containing: 1) aCL and/or aPS and $\alpha\beta_2$ GPI antibodies; or 2) moderate to strong levels (> 50 Units) of $\alpha\beta_2$ GPI in the absence of aCL or aPS antibodies were classified as APL positive; all other samples were considered APL negative. The (%) agreement was calculated by comparing the positive and negative results obtained by each assay (aPT, aCL, aPS, and $\alpha\beta_2$ GPI) with the Antiphospholipid (APL) status of the sample. For example, 100% agreement would mean that individual test results (positive or negative) were consistent with the APL status (positive or negative) in all samples of a given population.

% agreement with APL status

	IgM aPT	IgM aCL	IgM aPS	IgM anti- β_2GPI
Healthy - Serum	97.5%	92.5%	97.5%	97.5%
Autoimmune	80%	62.5%	72.5%	97.5%
Unselected SLE	82.5%	90%	80%	92.5%
Primary APS	81.8%	90.9%	72.7%	100%
Healthy - Plasma	97%	100%	99%	99%
LA + (group 1)	60%	75%	75%	90%
LA controls	87.5%	94.6%	98.2%	96.4%

Precision

Samples with known IgM aPT antibody levels (one low, one moderate, and one high) were assayed in 23 replicates on three different plates from two different kit lots. The mean intra-assay and inter-assay coefficients of variation (%CVs) are presented in the following table. The reported intra-assay coefficient of variation is the mean of six separate %CVs. The reported inter-assay %CV is the coefficient of variation obtained from the mean %CV of three plates from two lots.

IgM aPT Value Range	Mean Intra-assay %CV	Mean Inter-assay %CV
Low (< 5 M units)	7.2%	3.8%
Moderate (15-20 M units)	4.4%	1.5%
High (42-49 M units)	4.1%	2.5%

LIMITATIONS OF THE TEST

Anti-prothrombin (aPT) antibodies are frequently found in the serum or plasma of Lupus Anticoagulant (LA) positive patients, however, not all aPT are associated with LA activity. The REAADS IgM anti-Prothrombin ELISA may detect aPT antibodies associated as well as those not associated with LA activity. The anti-prothrombin (aPT) antibody concentration values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures. If clinical findings suggest the presence of antiphospholipid antibodies and the patient is negative for antiprothrombin antibodies, some investigators recommend testing for anti-cardiolipin antibodies, anti-phosphatidylserine antibodies, anti- β_2 glycoprotein I antibodies, and the lupus anticoagulant (LA) to confirm the negative result. A patient may be considered positive for antiphospholipid antibodies if one or all of the tests give positive results.

Warranty

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone +1-303-457-4345, fax +1-303-457-4519, email techsupport@Corgenix.com, or contact a Corgenix authorized distributor.

**REAADS®
IgM Anti-Prothrombin Semi-Quantitative Test Kit****In-vitro-Diagnostikum**

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REFERENCES

1. Harris EN, Gharavi AE, Hughes GRV. Anti-phospholipid Antibodies. *Clin Rheum Dis.* 11:591, 1985.
2. Galli M, Comfurius P, et al. Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein cofactor. *Lancet* 335:1544, 1990.
3. McNeil HP, Simpson RJ, et al. Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: B₂- glycoprotein I. *Proc Nat Acad Sci USA* 87:4120, 1990.
4. Bevers EM, Galli. M, Barbui T, et al. Lupus anticoagulant IgM's (LA) are not directed to phospholipids only, but to a complex of lipid-bound human prothrombin. *Thromb Haemost* 66:629, 1991
5. Roubey RA. Autoantibodies to phospholipid-binding proteins: a new view of lupus anticoagulants and other "antiphospholipid" antibodies. *Blood* 84:2854, 1994.
6. Permpikul P, Rao LVM, Rapaport SI. Functional and binding studies of the roles of prothrombin and B₂ glycoprotein I in the expression of lupus anticoagulant activity. *Blood* 83:2878, 1994.
7. Hunt JE, Krilis SA. The fifth domain of B₂-glycoprotein I contains a phospholipid binding site (Cys281-Cys288) and a region recognized by anticardiolipin antibodies. *J Immunol.* 152:653, 1994.
8. Matsuura E, Igarashi Y, Fujimoto M, et al. Heterogeneity of anticardiolipin antibodies defined by the anticardiolipin cofactor. *J Immunol.* 148:3885, 1992.
9. Guerin J, Feighery C, et al. Antibodies to B₂-glycoprotein I--a specific marker for the antiphospholipid syndrome. *Clin Exp Immunol.* 109:304, 1997.
10. Dier K, Keedy K, Fink C, Lopez L. Correlation of IgM, IgM and IgA anti-B₂ glycoprotein I antibodies with thrombosis in SLE patients. *Lupus* 7:S185, 1998.
11. Gomez-Pacheco L, Villa AR, et al. Serum anti-B₂-glycoprotein-I and anticardiolipin antibodies during thrombosis in systemic lupus erythematosus patients. *Am J Med.* 106:417, 1999.
12. Arvieux J, Darnige, L, et al. Development of an ELISA for autoantibodies to prothrombin showing their prevalence in patients with lupus anticoagulants. *Thromb Haemost* 74:1120, 1995.
13. Asherson RA, Khamashta MA, et al. The "primary" antiphospholipid syndrome: Major clinical and serological features. *Medicine* 68:366, 1989.
14. Wilson WA, Gharavi AE, Koike T, et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. Report of an international workshop. *Arthritis Rheum.* 42:1309, 1999.
15. Horbach DA, Oort EV, et al. Lupus anticoagulant is the strongest risk factor for both venous and arterial thrombosis in patients with systemic lupus erythematosus. *Thromb Haemost.* 76:916, 1996.
16. Galli M, Barbui T, et al. Prothrombin as cofactor for antiphospholipids. *Lupus* 7:suppl 2, S37, 1998.
17. Horbach DA, van Ort E, et al. The contribution of antiprothrombin-antibodies to lupus anticoagulant activity. *Thromb Haemost* 79:790, 1998.
18. Galli M, Barbui T, et al. Antiprothrombin antibodies: detection and clinical significance in the antiphospholipid syndrome. *Blood* 93:2149, 1999.
19. Galli M, Ruggeri L, Barbui T. Differential effects of anti-B₂-glycoprotein I and antiprothrombin antibodies on the anticoagulant activity of activated Protein C. *Blood* 91:1999, 1998.
20. Galli M, Finazzi G, Norbis F, et al. The risk of thrombosis in patients with lupus anticoagulants is predicted by their specific coagulation profile. *Thromb Haemost.* 81:695, 1999.
21. de Groot PG, Horbach DA, et al. Anti-prothrombin antibodies and their relation with thrombosis and lupus anticoagulant. *Lupus* 7:suppl 2, S32, 1998.
22. Guerin J, Smith O, White B, et al. Antibodies to prothrombin in antiphospholipid syndrome and inflammatory disorders. *Br J Haematol* 102:896, 1998.
23. Vaarala O, Puurunen M, et al. Antibodies to prothrombin imply a risk of myocardial infarction in middle aged men. *Thromb Haemost* 75:456, 1996.
24. Brandt JT, Triplett DA, et al. Criteria for the diagnosis of lupus anticoagulants: an update. *Thromb Haemost* 74:1185, 1995.
25. Data on File

SYMBOL LEGEND

										
Manufacturer	Authorized Representative	In vitro diagnostic medical device	Batch Code	Use by/ Expiry Date	Temperature Limitation	Irritant	Biological Risk	Catalog Number	European Conformity	Consult Instructions for Use/ Package Insert
Hersteller	Bevoll-mächtiger	In-vitro-Diagnostikum	Chargennummer	Verfallsdatum	Temperatur-beschränkungen	Reizend	Biologisches Risiko	Katalognummer	CE-Konformitäts-kennzeichnung	Gebrauchsanweisung im Inneren der Verpackung beachten
Fabriqué par	Représentant agréé	Dispositif de diagnostic in vitro	Code de Lot	Utiliser jusqu' à/ Date de péremption	Limites de température	Irritant	Risque biologique	Numéro de catalogue	Conformité aux normes européennes	Consulter le mode d'emploi/ notice jointe au conditionnement
Fabricado por	Representante autorizado	Dispositivo médico para diagnóstico in vitro	Código de Lote	Usar antes de/ Fecha de caducidad	Limitación de temperatura	Irritante	Riesgo biológico	Número de catálogo	Conformidad europea	Consultar las instrucciones de uso/ prospecto del envase
Prodotta da	Rappresentante autorizzato	Dispositivo medico-diagnostico in vitro	Codice del lotto	Scade il/ data di scadenza	Limite di temperatura	Irritante	Rischio biologico	Numero di catalogo	Conformità europea	Consultare le istruzioni per l'uso/ il foglietto illustrativo

MT Promedt Consulting GmbH
 Altenhofstrasse 80
 D-66386 St. Ingbert/Germany

Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020, USA

REAADS® is a registered trademark of Corgenix, Inc.
 © 2008, Corgenix, Inc.

10294 08
 Effective: 2008-01-11