**REAADS[®] VON WILLEBRAND
FACTOR ACTIVITY TEST KIT****For *In Vitro* Diagnostic Use****INTENDED USE**

An enzyme-linked immunosorbent assay (ELISA) for the quantitative determination of Von Willebrand Factor Activity (VWF:Act) in citrated human plasma.

SUMMARY AND EXPLANATION OF THE TEST

Von Willebrand Factor Antigen (VWF:Ag or Factor VIII-related protein) is a plasma protein found in circulation combined by non-covalent interactions with Factor VIII (FVIII:C), a pro-coagulant protein also known as the anti-hemophilic factor. These two proteins show distinct biochemical and functional properties as well as different antigenic determinants; their plasma levels may vary independently of each other.^{1,2} Deficiency of FVIII causes classic hemophilia while deficiency of VWF causes vonWillebrand disease. Most of VWF:Ag is synthesized and stored by endothelial cells while 15-20% is synthesized by megakaryocytes and stored in circulating platelets. A VWF:Ag unit has a molecular weight of about 250 kD and tends to polymerize in circulation, with multimers ranging in size from 850 kD to as large as 15×10^6 D.³

VWF plays a very important role in hemostasis; it protects FVIII from proteolytic cleavage in circulation and helps platelets to aggregate or to adhere to sites of vascular damage. The *in vivo* half-life of FVIII:C without VWF is shortened from 10-12 hours to a few minutes. These two mechanisms prevent bleeding. Von Willebrand disease is characterized by a deficiency or defect of VWF. Decreased VWF activity in plasma can be the result of low concentrations (quantitative or type I defect) or functional changes of VWF (qualitative or type II defect).^{4,5} Von Willebrand disease is the most common inherited bleeding disorder and is characterized by easy bruising and prolonged bleeding from mucosal surfaces. The prevalence of Von Willebrand disease has been estimated to be 1-3% of the general population. Greater than 70% of Von Willebrand disease patients have a type I deficiency while approximately 20% have a type II deficiency.⁶

The laboratory diagnosis of Von Willebrand disease may require both quantitative and qualitative (functional) determinations to differentiate the two predominant subtypes of the disease, type I and type II.^{7,8} The classification of Von Willebrand disease into subtypes is important in determining the course of clinical treatment.

PRINCIPLE OF THE TEST

The REAADS VWF:Act assay is a sandwich ELISA. A monoclonal capture antibody specific for the portion of VWF which binds platelets is coated to 96-microwell polystyrene plates.^{9,10} Diluted patient plasma is incubated in the wells, allowing any available antigen to bind to the monoclonal antibody on the microwell surface. The plates are washed to remove unbound proteins and other plasma molecules. Bound antigen is quantitated using horseradish peroxidase (HRP) conjugated anti-human VWF detection antibody. Following incubation, unbound conjugate is removed by washing. A chromogenic substrate of tetramethylbenzidine (TMB) and hydrogen peroxide (H_2O_2) is added to develop a colored reaction. The intensity of the color is measured in optical density (O.D.) units with a spectrophotometer at 450 nm. Patient VWF:Act in relative percent concentration is determined against a curve made from the reference plasma provided with the kit.

REAGENTS

Store at 2-8°C. Do Not Freeze.

Each REAADS Von Willebrand Factor Activity (VWF) 96-microwell Test Kit contains the following reagents:

- 12 x 8 stabilized antibody coated microwells (12 strips of 8 breakaway wells), with frame. Wells are coated with a monoclonal antibody specific for the functional region of VWF. (mouse)
- 60 mL Sample Diluent (blue-green solution); *
- 3 x 0.5 mL lyophilized Reference Plasma for preparation of reference curve, with assay sheet.
- 12 mL HRP-Conjugated anti-human VWF Antibody Solution (blue solution). (leporine)
- 13 mL Substrate (TMB/H₂O₂).
- 15 mL Stopping Solution (0.36N sulfuric acid).
- 30 mL Wash Concentrate [33X phosphate buffered saline (PBS) with 0.01% Tween 20]. Note: turbidity may appear in wash concentrate which will not affect component performance and should disappear when working dilution is prepared.

* **CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use

1. Human source material used to prepare the reference plasma included in this kit has been tested and shown to be negative for antibodies to HBsAg, HCV and HIV-I & II by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Substrate can cause irritation to the eyes and skin. Absorption through the skin is possible. Use gloves when handling substrate and wash thoroughly after handling. Keep reagent away from ignition sources. Avoid contact with oxidizing agents.
6. Certain components are labeled with the following:

Irritating to eyes (R 36). Avoid contact with skin (S 24). Avoid contact with eyes (S 25). In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S 26). If swallowed, seek medical advice immediately and show this container or label (S 46).

Irritant . Biological Risk .

SPECIMEN COLLECTION AND PREPARATION

Plasma collected with either 3.2% or 3.8% sodium citrate as an anticoagulant should be used as the sample matrix. Blood should be collected by venipuncture, and the sample centrifuged immediately. Remove the plasma and store at 2-8°C until testing can be performed. If not tested within 8 hours of collection, the sample should be stored at -70°C and tested within 1 month.

INSTRUCTIONS FOR USE

Materials Provided:

REAADS Von Willebrand Factor Activity Test Kit; see "Reagents," for a complete list.

Materials Required but not Supplied:

- VWF:Act Control Plasma. Reconstitute Control Plasma selected for use following manufacturer's instructions, and store as recommended.
- Reagent grade water (1 L) to prepare PBS/Tween 20 wash solution, to reconstitute Reference Plasma, and to zero or blank the plate reader during the final assay step.
- Graduated cylinders
- Precision pipettors capable of delivering between 5 and 1000 μL , with appropriate tips
- Miscellaneous glassware appropriate for small volume handling
- Flask or bottle, 1 liter
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semi-automated washing system
- Disposable gloves, powder-free recommended
- Plate reading spectrophotometer capable of reading absorbance at 450 nm (with a 650 nm reference, if available)
- Multichannel pipettors capable of delivering to 8 wells simultaneously
- Microdilution tubes for patient sample preparation

Procedural Notes

1. Bring plasma samples and kit reagents to room temperature (18–26°C) and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage (2–8°C) as soon as possible.
2. All dilutions of reference plasma, control plasma selected for use, and patient samples must be made just prior to use in the assay.
3. A single water blank well should be set up on each plate with each run. No sample or kit reagents are to be added to this well. Instead, add 200 μL of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to zero or blank against this water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution from a plastic squeeze bottle with a wide tip into the bottom of the microwells. Wash solution in the water blank well will not interfere with the procedure. An automated microtiter plate washing system can also be used.
5. **IMPORTANT:** Failure to adequately remove residual PBS/Tween 20 can cause inconsistent color development of the substrate solution.
6. Use a multichannel pipettor capable of delivering to 8 wells simultaneously when possible. This speeds the process and allows for more uniform incubation and reaction times for all wells.
7. Carefully controlled timing of all steps is critical. All reference plasma dilutions, controls and samples must be added to the microwells within a five minute period. Batch size of samples should not be larger than the amount that can be added within this time period.
8. For all incubations, the start of the incubation period begins upon the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below normal room temperature (18–26°C) may contribute to inaccurate results.
11. Avoid contamination of reagents when opening and removing aliquots from the primary vials.
12. Do not use kit components beyond expiration date.
13. Coated microwells, conjugate, and substrate are lot specific components that should not be used with different kit lots.

REAGENT PREPARATION

1. **Wash Solution (PBS/Tween 20):** Measure 30 mL Wash Concentrate (33X PBS/Tween 20) and dilute to 1 L with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused PBS solution in the refrigerator at 2–8°C. Discard if the solution shows signs of microbial or cross-contamination.
2. Reconstitute Reference Plasma by adding 0.5 mL reagent grade water. Swirl gently to mix. Allow to stand for 10 minutes before use for complete dissolution. Stable for 8 hours when stored at 2–8°C. Reconstitute appropriate control plasma following manufacturer's instructions, and store as recommended.

Assay Procedure

1. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
2. Assay each reference plasma dilution in duplicate. Duplicate determinations are also recommended for patient and control samples. One well should be run as a reagent blank; sample diluent without plasma is added to the well as explained in step 6 of this section. This well is treated the same as a patient sample in subsequent assay steps. A water blank well should be included with each plate; it is to remain empty until 200 µL of reagent grade water is added at the completion of the assay, immediately prior to reading the plate. The water blank well is to be used to zero the plate reader.
3. Using the Reference Plasma provided with the kit, prepare five reference plasma dilutions as described below.

Volume Reference Plasma		Volume Sample Diluent	*Reference Level (%)
100 µL	+	1000 µL =	200
50 µL	+	1000 µL =	100
25 µL	+	1000 µL =	50
25 µL	+	2000 µL =	25
25 µL	+	4000 µL =	12.5

* Reference level value to be used for constructing reference curve only.

4. Prepare a 1:21 dilution of each patient sample and control plasma selected for use in Sample Diluent (blue-green solution); e.g. 25 µL sample added to 500 µL Sample Diluent. Mix thoroughly.
5. Add 100 µL of the dilutions (reference plasmas x 5, patient samples and controls) to the appropriate microwells.
6. Add 100 µL of Sample Diluent to the reagent blank well. Leave the water blank well empty.
7. Incubate 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and empty the fluid. Do not allow samples to contaminate other microwells.
8. Wash 4 times with working wash solution (PBS/Tween 20). Each well should be filled with wash solution per wash. Wash solution in the empty well intended to serve as a water blank will not interfere with the procedure. Invert microwells between each wash to empty fluid. Use a snapping motion of the wrist to shake the liquid from the wells. The frame must be squeezed at the center on the top and bottom to retain microwell modules during washing. Blot on absorbent paper to remove residual wash fluid. Do not allow wells to dry out between steps.
9. Add 100 µL HRP-Conjugated Antibody Solution (blue) to each well (except the water blank well).
10. Incubate for 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and empty the conjugate solution.
11. Wash 4 times with working wash solution (PBS/Tween 20) as in step 8. Wash solution in the water blank well will not interfere with the procedure. Use a snapping motion to drain the liquid, and blot on absorbent paper after the final wash. Do not allow the wells to dry out.

12. Add 100 μ L Substrate to each well (except for the water blank well) and incubate for 10 minutes at room temperature. Add the substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
13. Add 100 μ L Stopping Solution (0.36N sulfuric acid) to each well (except for the water blank well) to stop the enzyme reaction. Be sure to add Stopping Solution to the wells in the same order and at the same rate as the Substrate was added. Blue substrate will turn yellow and colorless substrate will remain colorless. Do not add Stopping Solution to the water blank well. Instead, add 200 μ L reagent grade water to the water blank well. Blank or zero the plate reader against the water blank well. Read the O.D. of each well at 450 nm, against a 650 reference filter (if available). For best results, the O.D. values should be measured within 30 minutes after the addition of Stopping Solution.

RESULTS

Plasma VWF activity values are expressed as relative percent (%) as compared to pooled normal plasma.

1. Calculate the mean O.D. values for the duplicates of the reference plasma dilutions, controls selected for use, and patient samples.
2. Plot the mean O.D. obtained for each dilution of the reference plasma (x axis) against the corresponding value of the reference level (y axis). It is recommended that the curve be plotted on a log-log graph.
3. Using the mean O.D., determine the control and patient relative values from the graph, or alternatively, use linear regression to calculate from the reference curve.
4. To calculate VWF:Act levels in percent (%) of normal, multiply the control and patient relative values (obtained from the reference curve) by the assigned value for the REAADS Reference Plasma (see vial label).

For example: Patient relative value (from the reference curve): 40
Reference Plasma assigned value (from vial label): 105% of normal
Actual patient VWF:Act value (as % of normal): $40 \times 1.05 = 42\%$

5. Ensure that all quality control parameters have been met (see Quality Control) before reporting test results.

QUALITY CONTROL

1. The mean O.D. of the reagent blank should be less than 0.1 when the spectrophotometer has been blanked against the water well. Readings greater than 0.1 may indicate possible reagent contamination or inadequate plate washing.
2. O.D. values for the duplicates of the controls or patient samples should be within 20% of the mean O.D. value for samples with absorbance readings greater than 0.200.
3. VWF:Act values obtained for the controls should fall within manufacturer's assigned ELISA ranges. Occasional small deviations outside these ranges may be acceptable.
4. Each laboratory should periodically determine their own reference range for this assay.

EXPECTED VALUES¹³

One hundred normal plasma samples were tested on three lots of the REAADS VWF:Act assay. The mean value was 99.3% (SD of 39.0%) with 90% of the values for this normal population >50%. The assay cut-off was established at 50%. The 100 normal plasma samples were obtained from healthy individuals of unknown blood type. Individuals with type "O" blood type have been shown to have lower plasma levels of VWF (approximately 25% lower) when compared to those with other blood types. It is recommended that users of the REAADS Von Willebrand Factor Activity Kit establish reference ranges for the population served by their laboratory.

PERFORMANCE CHARACTERISTICS¹³

Detection Range:

The detection range for the REAADS VWF:Act assay has been determined to be 3.1-245.8%. This range is consistent with that published in the literature and reported by other commercially available assays.⁴⁻⁷ For greatest accuracy, samples that generate absorbance readings outside the O.D. range of the reference curve should be retested at an appropriate dilution. The lower limit of detection was calculated using data for the reagent blank plus 3 standard deviations, using 3 kit lots. The effective range of each run will depend on the assayed value of the reference plasma.

Precision:

Intra-assay precision:

Variability within a plate (intra-assay precision) was determined by testing three commercially prepared plasma samples across three plates from each of three lots (16 replicates per sample per plate). The data, presented in the following table, shows a mean CV of 4.3%. Intra-assay CVs were also determined by testing 14 commercially prepared plasma samples, in duplicate dilutions, across three plate lots. The mean CV for duplicates over three plate lots was 5.2%.

Inter-assay precision:

The precision between lots was determined by comparing the values recovered for eight different control samples on two lots. Each of the eight samples was tested in duplicate in each assay. The mean inter-assay CV was 8.3%.

Intra-assay precision (variability within a plate)	Von Willebrand Factor Activity range	CV range	Overall mean CV:
Replicates (x16) (3 lots):	(% of normal)		
	108.1-117.5%	2.7-5.0%	4.3%
	58.1-62.1%	2.6-6.4%	
37.3-38.4%	2.9-6.4%		
Duplicates (3 lots):	3.2-122.4%	0.7-16.5%	5.2%
Inter-assay inter-lot precision (variability between lots)			
Duplicates (2 lots):	4.8-116.4%	0.0-28.4%	8.3%

Curve Fit:

Curve Fit was evaluated by analyzing reference curves prepared from the Reference Plasma supplied with each kit. Reference Plasma dilutions were prepared and tested in duplicate (as outlined in the product insert) on two different plates, on three lots. For each assay, the mean O.D. of each Reference Plasma dilution was plotted against its Von Willebrand Factor activity concentration, in log-log format (log O.D. vs log of % Von Willebrand Factor activity concentration). The fit of each curve was calculated by linear regression, and expressed as the coefficient of determination (R-squared). The average coefficient of determination (R-squared) was 0.99296.

Percent Recovery:

Percent Recovery was determined by mixing control plasmas, with known values, in various proportions to assess how accurately the assay would recover theoretical values. The mean percent recovery was 102.1% based on two dilutions for each sample across three pilot lots, with an average percent variation of 3.4%.

LIMITATIONS OF THE TEST

The VWF:Act levels obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures. There is a normal plasma fluctuation of VWF due to unknown mechanisms. For this reason, repeat testing may be necessary. In addition, VWF acts as an acute phase reactant; it may be increased in various stressful conditions and diseases including pregnancy, oral contraceptive use, surgery, liver and autoimmune diseases, prostate cancer, etc.^{4,5}

Plasma samples can be inadvertently depleted or degraded of VWF by improper collection or laboratory processing. Individuals with "O" blood type have been shown to have lower plasma levels of VWF (25%) when compared to those with other blood types. Acquired Von Willebrand disease has been reported in some patients with lymphoproliferative disease.⁷

As with any assay employing antibodies from an animal source (e.g. mouse, rabbit, goat, etc.) to capture a target molecule, the possibility exists for interference in the serum or plasma of patients who have been exposed to preparations containing animal antibodies for diagnosis or therapy. Falsely elevated or depressed values may be seen in these patients.

Testing samples containing hemoglobin, lipids and bilirubin is not recommended as they may interfere with VWF Activity results.

The presence of Rheumatoid Factor in plasma samples can sometimes interfere in sandwich ELISAs by binding to the capture and/or detecting antibodies. Rheumatoid Factor interference has not been evaluated in the VWF Activity ELISA and should be considered when evaluating results.

Warranty

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone (303) 457-4345, fax (303) 457-4519, e-mail: technicalsupport@corgenix.com, or contact a Corgenix authorized distributor.

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REFERENCES

- Hoyer LW. The Factor VIII Complex: Structure and Function. *Blood* 58:1-12,1981.
- Lollar P. The Association of Factor VIII with Von Willebrand Factor. *Mayo Clinic Proceedings* 66:524-534, 1991.
- Ingerslev J. Von Willebrand Factor, Factor VIII and the Factor VIII/Von Willebrand Factor Complex. *Danish Medical Bulletin* 37:385-397, 1990.
- Ruggeri ZM. Structure and Function of Von Willebrand Factor: Relationship to von Willebrand's Disease. *Mayo Clinic Proceedings* 66:847-861, 1991.
- Montgomery RR, Collier BS. Von Willebrand Disease. In *Hemostasis and Thrombosis: Basic Principles and Clinical Practice*, p 134-163, 3rd Edition, JB Lippincott, Philadelphia, 1994.
- Ginsburg, D. Molecular Genetics of Von Willebrand Disease. *Thrombosis and Haemostasis* 82(2):585-591, 1999.
- Triplet DA. Laboratory Diagnosis of von Willebrand's Disease. *Mayo Clinic Proceedings* 66:832-840, 1991.
- Konkle BA. Laboratory Evaluation of Von Willebrand Disease. *Clin Chem* 41:489-490, 1995.
- Goodall AH. An immunoradiometric assay for human factor VIII/Von Willebrand Factor (VIII:VWF) Using a Monoclonal Antibody that Defines a Functional Epitope. *British Journal of Haematology* 59:565-577, 1985.
- Murdock PJ. Von Willebrand Factor Activity Detected in a Monoclonal Antibody-Based ELISA: an Alternative to the Ristocetin Cofactor Platelet Agglutination Assay for Diagnostic Use. *Thrombosis and Haemostasis* 78(4):1272-1277, 1997.
- Sadler, JE. Impact, Diagnosis and Treatment of Von Willebrand Disease. *Thrombosis and Haemostasis* 84(2):160-174, 2000.
- Gill JC, The effect of ABO Blood Group on the Diagnosis of von Willebrand Disease. *Blood* 69:1691-1695, 1987.
- Data on file.

SYMBOL LEGEND

										
Manufacturer	Authorized Representative	In vitro diagnostic medical device	Batch Code	Use by/ Expiry Date	Temperature Limitation	Irritant	Biological Risk	Catalog Number	European Conformity	Consult Instructions for Use/ Package Insert
Hersteller	Bevoll-mächtiger	In-vitro-Diagnostikum	Chargennummer	Verfallsdatum	Temperatur-beschränkungen	Reizend	Biologisches Risiko	Katalognummer	CE-Konformitäts-kennzeichnung	Gebrauchsanweisung im Inneren der Verpackung beachten
Fabriqué par	Représentant agréé	Dispositif de diagnostic in vitro	Code de Lot	Utiliser jusqu' à/ Date de péremption	Limites de température	Irritant	Risque biologique	Numéro de catalogue	Conformité aux normes européennes	Consulter le mode d'emploi/ notice jointe au conditionnement
Fabricado por	Representante autorizado	Dispositivo médico para diagnóstico in vitro	Código de Lote	Usar antes de/ Fecha de caducidad	Limitación de temperatura	Irritante	Riesgo biológico	Número de catálogo	Conformidad europea	Consultar las instrucciones de uso/ prospecto del envase
Prodotta da	Rappresentante autorizzato	Dispositivo medico-diagnostico in vitro	Codice del lotto	Scade il/ data di scadenza	Limite di temperatura	Irritante	Rischio biologico	Numero di catalogo	Conformità europea	Consultare le istruzioni per l'uso/ il foglietto illustrativo

MT Promedt Consulting GmbH
 Althenhofstrasse 80
 D-66386 St. Ingbert
 Germany

Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020 USA

REAADS® is a registered trademark of Corgenix, Inc.
 © 2006, Corgenix, Inc.

10755-03
 Effective: 2006-12-15