# READS ANA TEST

Cat. No. 10876:96 wells  
Cat. No. 10876:96 puits  
Art.-Nr. 10876:96 Auftragsstellen  
Cat. No. 10876:96 pocillos  
Cat. n. 10876:96 pozzetti

---

Corgenix, Inc.

11575 Main Street, Suite 400, Broomfield, Colorado 80020, USA  
Tel: 303-457-4345 Fax : 303-457-4519 [www.corgenixonline.com](http://www.corgenixonline.com)

- CONTENTS -  
- CONTENU -  
- INHALTSVERZEICHNIS -  
- CONTENIDO -  
- CONTENUTI -

English.....	3
Français.....	8
Deutsch.....	13
Español.....	19
Italiano.....	25
Symbols/Symboles/Symbole/Símbolos	
/Simboli.....	31

## ENGLISH

### INTENDED USE

The REAADS ANA TEST is semi-quantitative enzyme-linked immunosorbent assay (ELISA) for the detection of disease specific anti-Nuclear antibodies in human serum. The REAADS ANA TEST is intended for in vitro diagnostic use as an aid in the determination of autoimmune diseases.

### SUMMARY AND EXPLANATION

Anti-nuclear antibodies (ANA) are autoantibodies against various cell nucleus antigens and some of them are considered to be useful for diagnosis for autoimmune diseases such as Systemic Lupus Erythematosus (SLE), Mixed Connective Tissue Disease (MCTD), Sjögren's Syndrome (SjS), Progressive Systemic Sclerosis (PSS). Indirect immunofluorescence (IF) has been used widely to detect anti-nuclear antibodies. As the corresponding antigen can be presumed by staining pattern, it can yield information on a variety of autoantibody specificities. On the other hand, it is inconvenient when many samples are examined and it sometimes lacks assay specificity due to difference among reagents, operators and laboratory equipment.

The REAADS ANA TEST is a ELISA format test to detect disease specific anti-nuclear antibodies which are frequently observed in patients with autoimmune diseases. This test makes it possible to have objective interpretation since the results are expressed as Unit value and examines many samples in a short period of time.

### PRINCIPLE

The REAADS ANA TEST measures disease specific anti-Nuclear antibodies, anti RNP, Sm, SS-A, SS-B, Scl-70, Jo-1, CENP-B, Ribosomal-P, DNA and Histone antibodies present in the serum by ELISA. Calibrators and patient serum are added to microwell coated with antigens, allowing anti- Nuclear antibodies to react with the immobilized antigens (Sample incubation). After wash to remove any unbound serum proteins, horseradish peroxidase conjugated anti human IgG, IgA and IgM are added and incubated (Conjugate incubation). Following another washing step, the peroxidase substrate is added and incubated for an additional period of time (Substrate incubation). Acid solution is then added to each well to terminate the enzyme reaction and to stabilize the color development. The assay can be quantified by measuring the reaction photometrically and plotting the results.

### BRIEF ASSAY PROCEDURE

<Sample incubation> (20-25 °C) 60 min.	Add 100 µl of diluted sample (1:101) to each well of microwell plate ↓ Wash ↓
<Conjugate incubation> (20-25 °C) 60 min.	Add 100 µl of conjugate solution to each well ↓ Wash ↓
<Substrate incubation> (20-25 °C) 30 min	Add 100 µl of substrate to each well ↓ Add 100 µl of stop solution to each well ↓ Read absorbance ↓ Interpretation of result

## REAGENTS AND STORAGE

### 1) MICROWELL STRIPS

96 wells MICROWELL STRIPS (12 x 8 wells) coated with antigen produced from a mixture of recombinant purified proteins of Jo-1, RNP, SS-A, SS-B, CENP-B, Ribosomal P, in vitro transcribed U1 RNA, native purified proteins of SS-A, Sm, Scl-70, DNA and Histone, the breakaway strips packaged in a strip holder and sealed in a foil envelope with desiccant, are stable at 2-8°C until the labeled expiration date.

### 2) CALIBRATOR 1 (0 U/ml)

One vial containing 1.5 ml of Assay Diluent including 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

### 3) CALIBRATOR 2 (100 U/ml)

One vial containing 1.5 ml of ANA positive human control serum with Assay Diluent including 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

### 4) CONJUGATE CONCENTRATE (101X)

One vial containing 300 µl of horseradish peroxidase conjugated goat anti-human IgG, IgA and IgM (heavy chain specific) at a 101X concentration. Stable at 2-8°C until labeled expiration date. Bottle provided to prepare conjugate working solution.

### 5) CONJUGATE DILUENT

One bottle containing 24 ml of HEPES, Proclin 150 and BSA. Ready to use. Stable at 2-8°C until labeled expiration date.

### 6) ASSAY DILUENT

Two 50 ml bottles containing MOPS, Tween 20, bovine serum and 0.1% sodium azide. Ready to use. Stable at 2-8°C until labeled expiration date.

### 7) WASH CONCENTRATE (10X)

One 100 ml bottle containing PBS and Tween 20 as a 10X concentrate. Stable at 2-8°C until labeled expiration date.

### 8) SUBSTRATE SOLUTION

One 20 ml bottle containing 3,3',5,5'-tetramethylbenzidine dihydrochloride/hydrogen peroxide (TMB/H<sub>2</sub>O<sub>2</sub>). Ready to use. Stable at 2-8°C until labeled expiration date.

### 9) STOP SOLUTION

One 20 ml bottle containing 1N Sulfuric acid. Ready to use. Stable at 2-8°C until labeled expiration date.

### 10) POSITIVE CONTROL SERUM

One vial containing 0.2 ml of ANA positive human control serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

### 11) NEGATIVE CONTROL SERUM

One vial containing 0.2 ml of ANA negative human control serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

## PRECAUTIONS

- (1) This product is for in vitro diagnostic use only.
- (2) Do not use kit components beyond the stated expiration dates.
- (3) Avoid contact of reagents with eyes, skin and clothing. Reagents on skin must be washed away with plenty of water. TMB contains irritant and Stop Solution consists of a 1N sulfuric acid, which is a poison and corrosive.

- (4) Calibrator 2, Positive control and Negative control are derived from human serum, in which HBs antigen, HCV antibody, HIV-1 and HIV-2 antibodies cannot be detected. No test method, however, can guarantee the absence of these or any other infectious agents. These reagents and all patient samples should be handled as if they are capable of transmitting AIDS, hepatitis or any other infectious diseases.
- (5) Calibrator 1, Calibrator 2, Positive control, Negative control and Assay Diluent contain sodium azide (0.1%) as a preservative and must be handled with caution - do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (6) Some kit components contain animal origin materials, which are from non-infectious animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (7) Matching lot numbers of Microwell strips, Conjugate and Calibrator 2 must be used together in the assay. Do not substitute reagents from other kits.
- (8) All reagents must be brought to room temperature (20 -25°C) before starting the assay.
- (9) Do not expose the kit to direct sun during assay and storage.
- (10) Avoid microbial and cross contamination of reagents or samples.
- (11) Incubation temperatures above or below normal room temperature (20-25°C), shorter or longer time periods of incubation and inaccurate dilution may give erroneous results.
- (12) The wells must be rinsed with Wash Solution properly enough to avoid false positive.
- (13) Carefully pipette not to foam each sample and reagent to avoid cross contamination between microwells.
- (14) All microwell strips, which are not immediately required, should be returned to the zip lock pouch, which must be carefully resealed to avoid moisture absorption.
- (15) Wash concentrate may become turbid at 2-8°C, which does not cause inconsistent results.
- (16) The ANA values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedure.

## **MATERIALS REQUIRED BUT NOT PROVIDED**

- Microplate reader (wavelength: 450 nm, 620 nm/reference)
- Multichannel micropipette (e.g. 100 µl - 300 µl) for dispensing conjugate, substrate, and stop solution.
- Single channel pipette (10 µl & 100 µl)
- Reagent reservoir
- Autowasher or wash bottle
- Deionized or distilled water
- One liter graduated cylinder for preparation of wash solution
- Test tubes for patient sample dilutions (e.g. 1000 µl)
- Disposable pipette tips
- Paper towels
- Basin and disinfectant
- Microplate cover

## **PROCEDURE**

### **■ PREPARATION OF REAGENTS**

- Bring all assay materials to room temperature (20-25 °C) prior to use.
- Microwell Strips: Remove required microwell strip from pouch and place them in the frame. Promptly return unused strips to refrigerated storage.
- Wash Solution: Prepare 1:10 dilution of The Wash Concentrate prior to use (ex. add 100 ml of Wash Concentrate to 900 ml of distilled water). The diluted wash solution is stable for 2 weeks at 2 -8°C.

- Conjugate Solution: The Conjugate Reagent must be diluted prior to use. Dilute the Conjugate Reagent 1:101 with Conjugate Diluent. The diluted Conjugate Solution must be used within a day. Use disposable reservoir to avoid contamination.
- Do not dilute Calibrator 1, Calibrator 2, Assay Diluent, Substrate and Stop Solution, which are ready-to-use.

## ■ PREPARATION OF SAMPLES

- Use fresh patient sera. If storage is needed, they should be frozen below -20°C for up to one month, below -70°C for longer storage. Do not repeat freezing and thawing.
  - \*In case stored below -20°C for more than 6 months or freezing and thawing repeatedly, nonspecific results are obtained because of IgG denaturing.
- Dilute each patient serum 1:101 by adding 10µl of serum to 1ml of Assay Diluent.
  - \*Diluted samples may be stored up to 3 days if refrigerated.
  - \*Assay Diluent may form precipitate, which does not cause inconsistent results.

## ■ ASSAY PROCEDURE

### STEP 1. (SAMPLE INCUBATION)

Using the multi-channel pipetor, transfer 100µl of each diluted sample, Positive and Negative Controls into the appropriate microwells of the antigen test plate. Add Calibrators directly to appropriate wells. (Do not dilute Calibrators.)

\* Incubation starts on pipetting to the antigen-coated microwells. Pipetting should be completed as quickly as possible.

Cover wells with a plate sealer and incubate for 60 minutes at room temperature (20-25 °C).

### STEP 2. (WASHING)

Aspirate or discard the well contents. Fill the well with Wash Solution and then completely aspirate or discard the contents. Wash 4 times. Tap the plate on a paper towel to remove any remaining Wash Solution. When autowasher is used, wash 4 times.

\* Each laboratory is recommended to confirm its own appropriate washing times and set- up.

\* Wash Solution should be used at 20-25°C.

### STEP 3. (CONJUGATE INCUBATION)

Pour Conjugated Reagent into a reservoir. Add 100 µl of the Conjugated Reagent to each well with multichannel pipette. Cover wells with the plate sealer and incubate for 60 minutes at room temperature (20-25°C).

### STEP 4. (WASHING)

Wash the microplate following the STEP 2 procedure.

### STEP 5. (SUBSTRATE INCUBATION)

Pour Substrate into a reservoir and pipette 100 µl of the Substrate to each well with multichannel pipette.

\* The reservoir should be different from the one, which was used for pouring conjugate solution. A new disposable reservoir should be used because Substrate is easily oxidized by metal ion.

\* The Substrate, once poured in a reservoir, should not be returned to the bottle.

Cover wells with the plate sealer and incubate for 30 minutes at room temperature (20-25°C).

### STEP 6. (STOP REACTION)

Pour Stop Solution into a reservoir. Pipette 100 µl of the solution to each well with multichannel pipette.

## ■ READING

Read the absorbance of each well at 450 nm. If a dual wave length platereader is available, set the test wavelength at 450 nm and the reference at 620 nm.

\* Reading should be done as quickly as possible after stopping the reaction.

\* Ensure that the bottom of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading the plate.

## ■ CALCULATION OF RESULT

$$\text{Unit value (U/ml)} = \frac{(A_{450}\langle\text{Sample}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)}{(A_{450}\langle\text{Calibrator 2}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)} \times 100$$

\* A<sub>450</sub> is abbreviation of absorbance value at 450 nm.

\* An appropriate international reference material for measuring anti nuclear antibody by ELISA is not available; the assay is calibrated in relative arbitrary units.

## ■ QUALITY CONTROL

Each assay result should meet the following criteria.

A<sub>450</sub> of Calibrator1: ≤0.100

A<sub>450</sub> of Calibrator 2: >0.500

The Positive and Negative Controls must give the following results:

	Positive Control	Negative Control
Anti-ANA value (U/ml)	>40	<21

If any of these are not met, the results are invalid and the test should be repeated.

Before repeating assay, check the following procedure.

- Incubation Temperature
- Incubation Period of Time
- Washing

## TEST INTERPRETATION AND EXPECTED VALUES

The following is intended only as a guide for interpretation. Each laboratory is recommended to establish its own criteria for test interpretation based on sample populations typically encountered.

Anti-ANA value (U/ml)	Interpretation
<21	Negative for anti-ANA Ab
Greater than or equal to 21	Positive for anti-ANA Ab

The above ANA value was determined by ROC analysis with 695 serum samples of autoimmune disease and healthy blood donors.

## ■ LIMITATIONS

As with other diagnostic test procedures, the results obtained with the REAADS TEST ANA serve only as an aid to diagnosis and should not be interpreted as diagnostics in themselves.

## PERFORMANCE CHARACTERISTICS

### ■ CLINICAL SPECIFICITY AND SENSITIVITY

Disease	Positive sample	Positive rate
MCTD	55/56	98.2%
SLE	96/102	94.1%
SjS	44/46	95.7%
SSc	52/55	94.5%
PM/DM	18/27	66.7%
Normal Serum	2/38	5.3%

MCTD: Mixed connective tissue disease

SLE: Systemic lupus erythematosus

SjS: Sjögren's Syndrome

SSc: Systemic sclerosis

## ■ PRECISION

Repeatability was demonstrated by testing 3 samples in sextuple. Reproducibility was determined by testing 3 samples on 5 different days (day-to-day) and by testing 3 samples in 3 lots (lot-to-lot). %CV values for reproducibility and repeatability were below 15% for each sample.

## ■ ASSAY RANGE

The assay range of this kit is from 5U/ml to 200U/ml.

## ■ INTERFERING SUBSTANCES

Hemoglobin (up to 470 mg/dl), Bilirubin (up to 17.7mg/dl) and/or chyle (up to 2,490 unit as Formazine) are not affective on the assay result, but avoid using highly hemolysed samples or highly lipemic samples.

## REFERENCES

1. Tsay, G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. *Arthritis & Rheum.* 30: 389-396, (1987).
2. Tojo, T. et al.: Diversity of autoantibodies and their clinical significance in systemic autoimmune diseases. *Medical immunology* 14: 761-766, (1987).
3. Sharp, G.C. et al.: Mixed connective tissue disease: an apparently distinct rheumatic disease syndrome associated with a specific antibody to an extractable nuclear antigen (ENA). *Amer. J. Med.* 52: 148-152, (1972).
4. Query, C.C. et al.: A common RNA recognition motif identified within a defined U1 RNA binding domain of the 70K U1 snRNP protein. *Cell* 57: 89-101, (1989).
5. Query, C.C. et al.: A human autoimmune protein associated with U1 RNA contains a region of Homology that is cross-reactive with retroviral p30gag antigen. *Cell* 51: 211-220, (1987).
6. Krapf, A.R., C.A. von Mühlen, F.E. Krapf, R.M. Nakamura, E. Tan, *Atlas of Immunofluorescent Autoantibodies*, Urban & Schwarzenberg, Baltimore, MD. , 1996.

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1. Tsay, G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. *Arthritis & Rheum.* 30: 389-396, (1987).
2. Tojo, T. et al.: Diversity of autoantibodies and their clinical significance in systemic autoimmune diseases. *Medical immunology* 14: 761-766, (1987).
3. Sharp, G.C. et al.: Mixed connective tissue disease: an apparently distinct rheumatic disease syndrome associated with a specific antibody to an extractable nuclear antigen (ENA). *Amer. J. Med.* 52: 148-152, (1972).
4. Query, C.C. et al.: A common RNA recognition motif identified within a defined U1 RNA binding domain of the 70K U1 snRNP protein. *Cell* 57: 89-101, (1989).
5. Query, C.C. et al.: A human autoimmune protein associated with U1 RNA contains a region of Homology that is cross-reactive with retroviral p30gag antigen. *Cell* 51: 211-220, (1987).
6. Krapf, A.R., C.A. von Mühlen, F.E. Krapf, R.M. Nakamura, E. Tan, *Atlas of Immunofluorescent Autoantibodies*, Urban & Schwarzenber, Baltimore, MD. , 1996.

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### BIBLIOGRAFIA

1. Tsay, G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. *Arthritis & Rheum.* 30: 389-396, (1987).
2. Tojo, T. et al.: Diversity of autoantibodies and their clinical significance in systemic autoimmune diseases. *Medical immunology* 14: 761-766, (1987).
3. Sharp, G.C. et al.: Mixed connective tissue disease: an apparently distinct rheumatic disease syndrome associated with a specific antibody to an extractable nuclear antigen (ENA). *Amer. J. Med.* 52: 148-152, (1972).
4. Query, C.C. et al.: A common RNA recognition motif identified within a defined U1 RNA binding domain of the 70K U1 snRNP protein. *Cell* 57: 89-101, (1989).
5. Query, C.C. et al.: A human autoimmune protein associated with U1 RNA contains a region of Homology that is cross-reactive with retroviral p30gag antigen. *Cell* 51: 211-220, (1987).
6. Krapf, A.R., C.A. von Mühlen, F.E. Krapf, R.M. Nakamura, E. Tan, *Atlas of Immunofluorescent*
7. *Autoantibodies*, Urban & Schwarzenberg, Baltimore, MD. , 1996.

\*\*The following symbols are used on the label.

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**Manufactured for:**

Corgenix, Inc.  
11575 Main Street, Suite 400  
Broomfield, Colorado 80020 USA

**REAADS**<sup>®</sup> is a registered trademark of Corgenix, Inc.

© 2005, Corgenix, Inc

10930 03  
Effective: 2006-09-27

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MT Promedt Consulting GmbH