Anti-Polymer Antibody (APA) ELISA

For *In Vitro* Diagnostic Use

An enzyme-linked immunosorbent assay (ELISA) for the semi-quantitative determination of IgG anti-polymer antibodies in human serum.

INTENDED USE

The patented APA Test Kit is a semi-quantitative ELISA which detects IgG anti-polymer antibodies (APA) and identifies certain fibromyalgia patients. The APA Test Kit is intended for use: 1) as an aid in the diagnosis of patients presenting with the symptoms and signs of fibromyalgia syndrome, 2) as an aid in differentiating fibromyalgia patients from patients with other autoimmune diseases, and 3) as an aid in identifying a subgroup of fibromyalgia patients having a unique characteristic immune response.

SUMMARY AND EXPLANATION

Fibromyalgia Syndrome (FMS) is a chronic disorder characterized by subjective musculoskeletal pain and semi-objective tenderness to deep palpitation at defined somatic tender points. The American College of Rheumatology (ACR) criteria for the classification of FMS have provided a highly sensitive and specific means of identifying affected individuals. Based on those criteria, estimates for the prevalence of FMS in the general population ranges from 2 to 10% with females being affected 7 – 9 times more commonly than males.

A number of systemic manifestations occurring as co-morbidity in some FMS patients, but not required by the ACR criteria, include chronic widespread pain, insomnia, fatigue, morning stiffness, weakness, irritable bowel, interstitial cystitis, and reactive depression. One important benefit derived from establishing reliable criteria for the diagnosis of FMS has been the relatively homogeneous group of patients that can be assembled for a systematic evaluation of proposed laboratory abnormalities.

A variety of clinical similarities between Silicone Breast Implant (SBI) recipients and Fibromyalgia patients led to investigations of the APA levels in Fibromyalgia patients. Wilson, et al found that the APA positive rate in patients with FMS (47%) was higher than that in patients with osteoarthritis [OA] (19%) or rheumatoid arthritis [RA] (8%). The same was true for comparisons with lupus erythematosus [SLE] (3%), and progressive systemic sclerosis [PSS] (3%). The prevalence of APA seroreactivity was also significantly higher in patients with severe FMS (61%) than in FMS patients with less severe symptoms. These results suggest that the APA assay might be an objective marker for both the diagnosis and the clinical prognosis of Fibromyalgia.

PRINCIPLE OF THE TEST

The test is performed as an indirect ELISA. Diluted serum samples, calibrator, and controls are incubated in polymer coated microwells. Incubation allows the anti-polymer antibodies present in the samples to react with the immobilized coated surface. After the removal of unbound proteins by washing, antibodies specific for human IgG labeled with horseradish peroxidase (HRP) are added forming complexes with the polymer bound antibodies. Following another washing step, the bound enzyme-antibody conjugate is assayed by the addition of a single solution containing tetramethylbenzidine (TMB) and hydrogen peroxide (H₂O₂) as the chromogenic substrate. Color develops in the wells at an intensity proportional to the concentration of anti-polymer antibodies.

Results are obtained by reading the optical density (OD or absorbance) of each well on a spectrophotometer. A calibrator serum is provided with the IgG anti-polymer antibody concentration expressed in Units. Dividing the concentration value of the calibrator sera by the OD of the calibrator provides a conversion factor. The OD reading of the other samples are multiplied by the conversion factor to obtain IgG anti-polymer antibody concentration in Units.

REAGENTS

Store at 2 – 8°C . Do Not Freeze.

Each APA 96-microwell ELISA Kit contains the following reagents

- 96 APA Antigen Coated Microwells (12 strips of 8 wells), with frame
- 2 bottles (60 mL) APA Sample Diluent (green solution)
- 1 vial (0.250 mL) APA Calibrator (human source) *
- 1 vial (0.250 mL) APA Positive Control (human source) *
- 1 vial (0.250 mL) APA Normal Control (human source)*
- 1 bottle (15 mL) APA HRP-Conjugated Antibody Solution (goat anti-Human IgG)
- 1 bottle (15 mL) One Component Substrate (TMB/H₂O₂); ready to use
- 1 bottle (15 mL) Stopping Solution (0.36 N sulfuric acid)
- 2 bottles (30 mL) Wash Concentrate (33X PBS/Tween 20)

*** CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use

1. Human source materials used to prepare the calibrator and controls included in this kit have been tested and shown to be negative for antibodies to HBsAg, HCV, Syphilis, and HIV 1 & 2 by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Certain components of this product contain sodium azide as a preservative. Sodium azide has been reported to form lead and copper azides when left in contact with these metals. These metal azides are explosive. Any solutions containing azide must be thoroughly flushed with copious amounts of water to prevent the build-up of explosive metal azides in the plumbing system
6. One component substrate can cause irritation to the eyes and skin. Absorption through the skin is possible. Use gloves when handling substrate and wash thoroughly after handling. Keep reagent away from ignition sources. Avoid contact with oxidizing agents.
7. Certain components are labeled with the following: Harmful if swallowed (R 22). Irritating to eyes and skin (R 36/38). Avoid contact with skin and eyes (S 24/25). In case of contact with eyes, flush affected areas with copious amounts of water and seek medical advice (S 26). Wear suitable protective clothing (S 36).

SPECIMEN COLLECTION AND PREPARATION

Serum is the preferred sample matrix. Blood should be collected by venipuncture, and the serum separated from the cells by centrifugation after clot formation. If not tested immediately, specimens should be stored at 2-8°C . If specimens are to be stored for more than 72 hours, they should be frozen at -20°C or below. Avoid repeated freezing and thawing. Do not use hemolyzed, icteric, or lipemic serum as these conditions may cause aberrant results. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

INSTRUCTIONS FOR USE

Materials Provided:

Anti-Polymer Antibody (APA) ELISA; see "Reagents" for a complete list.

Materials Required but not Supplied:

- Reagent grade water to prepare PBS/Tween 20 wash solution (1L)
- Graduated cylinders
- Precision pipettors capable of delivering between 10 and 1000 μL , with appropriate tips
- Miscellaneous glassware or dilution tubes appropriate for small volume handling
- Flask or bottle, 1 liter
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semi-automated washing system
- Disposable gloves
- Plate reading spectrophotometer capable of reading absorbance at 450 nm (with a 650 nm reference, if available)
- Multichannel pipettors capable of delivering 100 μL to eight wells simultaneously
- Reagent reservoirs to facilitate delivery of reagents by multichannel pipettors

Procedural Notes

1. Bring serum samples and kit reagents to ambient temperature (18-26°C) and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage as soon as possible.
2. All dilutions of calibrator, control, and test sera must be made just prior to use in the assay.
3. An Air Blank (zero instrument on air prior to final reading of wells) or a Water Blank (optional) must be used on each plate with each run. No sample or kit reagents are to be added to the Water Blank well. For a Water Blank, add 200 μl of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to zero or blank against air or the water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution from a plastic squeeze bottle with a wide tip into the bottom of the microwells. Wash solution in the water blank well will not interfere with the procedure. An automated microtiter plate washing system can also be used.
5. **IMPORTANT:** Failure to adequately remove residual PBS/Tween 20 can cause inconsistent color development of the Substrate Solution.
6. Use a multichannel pipettor capable of delivering to eight wells simultaneously when possible. This speeds the process and provides more uniform incubation and reaction times for all wells.
7. Carefully controlled timing of all steps is critical. All calibrators, controls, and samples must be added to the plate within a five-minute period. Batch size of samples should not be larger than the number that can be added within this time period.
8. For all incubations, the start of the incubation period begins with the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below the ambient temperature of 18-26°C may contribute to inaccurate results.
11. Avoid contamination of reagents when opening and removing aliquots from the primary vials. Pour or pipette only the volume of reagents needed. Do not return excess volumes of either HRP-Conjugated Antibody Solution or One Component Substrate Solution to the original bottles.
12. Do not use kit components beyond expiration date.
13. Do not use kit components from different kit lot numbers.

REAGENT PREPARATION

Wash Solution (PBS/Tween 20): Measure 30 mL of Wash Concentrate (33X PBS/Tween 20) and dilute to 1 liter with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused wash solution in the refrigerator at 2 - 8°C. Discard if the solution shows signs of microbial or other contamination.

Assay Procedure

1. Assay the calibrator and controls in duplicate. It is advised that duplicate determinations also be made for each patient sample. At least one well should be run as a reagent control; APA Sample Diluent without serum or plasma is added to a well. This well will be treated the same as sample wells in subsequent assay steps.
2. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
3. Prepare a 1:101 dilution of the calibrator, controls, and patient samples in APA Sample Diluent (green solution), e.g., 10 μ L sample plus 1000 μ L APA Sample Diluent equals a 1:101 sample dilution.
4. Add 100 μ L of diluted calibrator, controls, and patient samples to the appropriate microwells.
5. Add 100 μ L of APA Sample Diluent to the reagent control well.
6. Incubate 60 minutes at 18-26°C. After the incubation is complete, carefully invert the microwells and dump or aspirate the sample fluid. Do not allow samples to contaminate other microwells.
7. Wash 4 times with PBS/Tween 20 wash solution. Grasp the plate frame and squeeze at the top and bottom to retain microwells during washing. Use a snapping motion of the wrist to shake the liquid from the wells. Completely fill each well with wash solution, then remove or shake out the wash solution into a container or sink Repeat the process three additional times. Invert the wells on absorbent paper and tap to blot residual wash solution. Do not allow the wells to dry out at any time.
8. Add 100 μ L HRP-Conjugated Antibody Solution (blue solution) to all wells. Incubate for 30 minutes at 18-26°C. After the incubation is complete, carefully invert the microwells and dump or aspirate the conjugate solution.
9. Wash 4 times with PBS/Tween 20 wash solution as in step 7.
10. Add 100 μ L One Component Substrate Solution to each well and incubate for 30 minutes at 18-26°C. Add the substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
11. Add 100 μ L Stopping Solution (0.36 N sulfuric acid) to each well to stop the enzyme reaction. Be sure to add the Stopping Solution to the wells in the same order and at the same rate as the substrate was added. Blue substrate color reactions will turn yellow and colorless solution will remain colorless when Stopping Solution is added.
12. If a Water Blank is used (see Procedural Notes #3), place 200 μ L of reagent grade water to the well designated as the water blank. Blank or zero the plate reader against air or a water blank well. Read the OD of each well at 450 nm (and 650 nm reference, if dual beam). The OD values should be measured within 30 minutes after the addition of Stopping Solution.

Results

1. Calculate the mean OD values for the duplicates of Calibrator, Controls, and patient samples.
2. Divide the assigned value of the Calibrator (printed on the vial label) by the mean OD of the Calibrator serum to obtain the conversion factor.
3. Multiply the mean OD for each of the Controls and patient samples by the conversion factor to obtain the Units (U) for each sample.

$$\text{Conversion Factor} = \frac{\text{Assigned Calibrator Value}}{\text{Mean Calibrator OD}}$$

$$\text{Sample Unit (U) Value} = \text{Conversion Factor} \times \text{OD of the sample.}$$

4. The Conversion Factor must be calculated for each assay run. Using a Conversion Factor from another assay will invalidate the results.
5. If the sample OD is higher than the upper limit of the plate reader, the sample should be further diluted 1:10 in APA Sample Diluent (i.e., add 10 μ L sample to 90 μ L APA sample diluent). Then proceed with the 1:101 dilution and re-assay to obtain a value. Multiply the resulting assay Unit by 10 to obtain the calculated Unit value.

Quality Control

1. The mean OD value of the Calibrator should be at least 0.400 to assure that the kit is functioning properly. Calibrator OD readings of less than 0.400 may indicate that the kit is no longer suitable for use.
2. The mean OD of the reagent control should be less than 0.100. Readings greater than 0.100 may indicate possible reagent contamination or inadequate plate washing.
3. The APA Unit (U) obtained for the Control Sera should be within the ranges indicated on the vial labels.
4. OD values for the duplicates of the Controls or patient samples should be within 20% of each other for samples with OD/absorbance readings greater than 0.200. If not, re-run the patient sample in duplicate.

NORMAL RANGE

Serum samples from 100 healthy blood donors were tested for IgG anti-polymer antibodies. The normal range was established as less than 30 Units (U).

PERFORMANCE CHARACTERISTICS

Clinical Specificity

Random Donor Controls:

Serum samples from 100 random blood donors were assayed for the presence of anti-polymer antibodies. Using the established cut-off of 30 Units (U), this normal population showed a specificity of 75%. This rate correlates with a previously published study² where the presence of anti-polymer antibodies were found in 19% of random donor controls.

Clinical Sensitivity:

Systemic Lupus Erythmatosis (SLE):

Serum samples from 93 patients diagnosed with Systemic Lupus Erythmatosis (SLE) that were collected in sequential order at a Rheumatology clinic were analyzed in the APA ELISA to determine the prevalence of anti-polymer antibodies. The presence or absence of FMS for each patient was not determined. The mean value was 16.4 Units (U) with a positive rate of 11.8%.^{2,5}

Rheumatoid Arthritis (RA):

Serum samples from 71 patients diagnosed with Rheumatoid Arthritis (RA) that were collected in sequential order at a Rheumatology clinic were analyzed in the APA ELISA to determine the prevalence of anti-polymer antibodies. The presence or absence of FMS for each patient was not determined. The mean value was 31.9 Units (U) with a positive rate of 28.2%.^{2,6}

Fibromyalgia Syndrome (FMS):

Serum samples from 16 patients diagnosed with Fibromyalgia Syndrome (FMS) were analyzed in the APA ELISA to determine the prevalence of anti-polymer antibodies. The mean value was 54.3 Units (U) with a positive rate of 43.8%.²

Summary:

	Random Donors	SLE	RA	FMS
#of Samples (n)	100	93	71	16
Mean (Units)	36.9	16.4	31.9	54.3
% Positive	25%	11.8%	28.2%	43.8%

Precision:

Two samples with known IgG APA values (one high and one low) were assayed in 23 replicates on three different plates. The mean intraassay and interassay percent coefficients of variation (%CVs) are presented in the following table. The reported intraassay coefficient of variation is the mean of the three separate intraassay CVs. Interassay CV is the coefficient of variation obtained from three plates from one lot.

IgG APA Value	Mean Intra-assay CV	Mean Inter-assay CV
High (254 Units)	6.7%	6.6%
Low (40 Vnits)	7.5%	4.0%

LIMITATIONS OF THE TEST

The IgG APA concentration values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures.

WARRANTY

This product is warranted to perform as described in this package insert. Corgenix Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661, Outside the United States phone 1-303-457-4345, fax 1-303-457-4519, email techsupport@Corgenix.com, or contact a Corgenix authorized distributor.

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REFERENCES

1. Xiao YM, Russell IJ, Michalek JE, Wilson RB. Anti-Polymer Antibodies Identify a Large Subgroup of Fibromyalgia Syndrome Patients. Presentation; Myopain 2004.
2. Wilson RB, Gluck OS, Tesser JR, Rice JC, Meyer A, Bridges AJ. Anti-Polymer Antibody Reactivity in a Subset of Patients with Fibromyalgia Correlates With Severity. J Rheumatol 1999;26(2):402-407.
3. Bennett RM. Fibromyalgia and the Disability Dilemma. A New Era in Understanding a Complex, Multidimensional Pain Syndrome. Arthritis Rheum 1996;39(10):1627-1634.
4. Wolfe F, Smythe HA, Yunus MB, Bennett RM, Bombardier C, Goldenberg DL, et al. The American College of Rheumatology 1990 Criteria for the Classification of Fibromyalgia. Report of the Multicenter Criteria Committee. Arthritis Rheum 1990;33(2):160-172.
5. Gladman, D.D. Urowitz M.B. Gough, J., MacKinnon, A. Fibromyalgia is a Major Contributor to Quality of Life in Lupus. J Rheumatol, 1997. 24(11): p. 2145-8
6. Wolfe, F. and K. Michaud, Severe Rheumatoid Arthritis (RA), Worse Outcomes, Co-Morbid Illness, and Sociodemographic Disadvantage Characterize RA Patients with Fibromyalgia. J Rheumatol, 2004. 31(4): p. 695-700

									
Authorized Representative	In vitro diagnostic medical device	Batch	Expiry Date	Storage conditions	Harmful	Irritant	Biological Risk	Catalog Number	European Conformity
(Bevollmächtigter)	(In-vitro-Diagnostikum)	(Charge)	(Verfallsdatum)	(Lagerbedingungen)	(Gesundheitsschädlich)	(Reizend)	(Biologisches Risiko)	Katalognummer	CE-Konformitätskennzeichnung
(Représentant agréé)	(Dispositif de diagnostic in vitro)	(Lot)	(Date de péremption)	(Conditions de stockage)	(Nocif)	(Irritant)	(Risque biologique)	Numero di catalogo	Conformità europea
(Representante autorizado)	(Dispositivo médico para diagnóstico in vitro)	(Lote)	(Fecha de caducidad)	(Condiciones de almacenamiento)	(Perjudicial)	(Irritante)	(Riesgo biológico)	Numéro de catalogue	Conformité aux normes européennes
(Rappresentante autorizzato)	(Dispositivo medico-diagnostico in vitro)	(Lotto)	(Data di scadenza)	(Condizioni di conservazione)	(Nocivo)	(Irritante)	(Rischio biologico)	Número de catálogo	Conformidad europea

MT Promedt Consulting GmbH
 Altenhofstrasse 80
 D-66386 St. Ingbert/ Germany

Manufactured by:

Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020, USA
 © 2007

The APA Assay is covered by U.S. Patent No. 5,834,215
 European Patent No.0873518
 Australian Patent No. 710252
 and other patents and pending patents.

11679-E 01
 Effective: 2007-02-23