READS anti-RNP

ELISA Test Kit

Cat. No. 10869: 96 wells
Cat. No. 10869: 96 puits
Art.-Nr. 10869: 96 Auftragsstellen
Cat. No. 10869: 96 pocillos
Cat. n. 10869: 96 pozzetti

Corgenix, Inc.

11575 Main Street, Suite 400, Broomfield, Colorado 80020, USA

Tel: 303-457-4345 Fax : 303-457-4519 www.corgenixonline.com

- CONTENTS -
- CONTENU -
- INHALTSVERZEICHNIS -
- CONTENIDO -
- CONTENUTI -

English.....	1
Français.....	7
Deutsch.....	13
Español.....	19
Italiano.....	25

Symbols/Symboles/Symbole/Símbolos /Simboli...31

English

INTENDED USE

The REAADS anti-RNP ELISA Test Kit is semi-quantitative enzyme-linked immunosorbent assay (ELISA) for the detection of anti-RNP antibodies in human serum. The REAADS anti-RNP ELISA Test Kit is intended for in vitro diagnostic use as an aid in the determination of certain autoimmune diseases.

SUMMARY AND EXPLANATION

RNP antibodies are ENA (extractable nuclear antigen) antibodies, present with high frequency in the sera of patients with collagen diseases such as SLE and mixed connective tissue disease (MCTD). Mixed connective tissue disease is a clinical disease combining features of SLE, Progressive Systemic Sclerosis (PSS) and Polymyositis/Dermatomyositis (PM/DM). The presence of RNP antibodies in the sera of patients is essential for diagnosis of MCTD. RNP antibodies, when present alone at high levels, are diagnostic of MCTD. Lower levels of anti-RNP, in conjunction with other autoantibodies may be observed in PSS, Sjögren's Syndrome (SS) and Rheumatoid Arthritis (RA).

Among the methodologies available to detect anti-ENA's are Immunoblot, Double Immunodiffusion (DID), Immunoprecipitation (IPP) and Enzyme-linked Immunosorbent Assay (ELISA). The ELISA tests are specific, sensitive, and due to their objectivity and rapidity, suitable for testing samples from patients with suspected connective tissue diseases.

The REAADS anti-RNP ELISA Test Kit uses recombinant protein for the solid phase antigen. Therefore, the test shows specificity for detection of anti RNP antibodies present in high frequency in the sera of patients with certain connective tissue diseases.

PRINCIPLE

The REAADS anti-RNP ELISA Test Kit measures anti-RNP antibodies present in the serum by ELISA. Calibrators and patient serum are added to microwell coated with RNP antigens, allowing anti-RNP antibodies to react with the immobilized antigen (Sample incubation). After wash to remove any unbound serum proteins, horseradish peroxidase conjugated anti human IgG, IgA and IgM are added and incubated (Conjugate incubation). Following another washing step, the peroxidase substrate is added and incubated for an additional period of time (Substrate incubation). Acid solution is then added to each well to terminate the enzyme reaction and to stabilize the color development. The assay can be quantified by measuring the reaction photometrically and calculating the results.

BRIEF ASSAY PROCEDURE

<Sample incubation> (20-25°C) 60 min.	Add 100 µl of diluted sample (1:101) to each well of microwell plate ↓ Wash ↓
<Conjugate incubation> (20-25°C) 60 min.	Add 100 µl of conjugate solution to each well ↓ Wash ↓
<Substrate incubation> (20-25°C) 30 min	Add 100 µl of substrate to each well ↓ Add 100 µl of stop solution to each well ↓ Read absorbance ↓ Interpretation of result

REAGENTS AND STORAGE

1) RNP MICROWELL STRIPS

96 wells MICROWELL STRIPS (12 x 8 wells) coated with antigens produced from recombinant purified proteins, the breakaway strips packaged in a strip holder and sealed in a foil envelope with desiccant, are stable at 2-8°C until the labeled expiration date.

2) CALIBRATOR 1 (0U/ml)

Two vials containing 1.5 ml of Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

3) CALIBRATOR 2 (100U/ml)

Two vials containing 1.5 ml each of anti RNP antibody positive human serum with Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

4) CONJUGATE REAGENT

One vial containing 15 ml of horseradish peroxidase conjugated goat anti-human IgG, IgA and IgM (heavy chain specific), HEPES, Proclin 150 and BSA. Ready to use. Stable at 2-8°C until labeled expiration date.

5) ASSAY DILUENT*

Two 50 ml bottles containing PBS, Tween 20 and 0.1% sodium azide. Ready to use. Stable at 2-8°C until labeled expiration date.

6) WASH CONCENTRATE (10X) *

One 100 ml bottle containing PBS and Tween 20 as a 10X concentrate. Stable at 2-8°C until labeled expiration date.

7) SUBSTRATE SOLUTION*

One 20 ml bottle containing 3,3',5,5'-tetramethylbenzidine dihydrochloride/hydrogen peroxide (TMB/H₂O₂). Ready to use. Stable at 2-8°C until labeled expiration date

8) STOP SOLUTION*

One 20 ml bottle containing 1N Sulfuric acid. Ready to use. Stable at 2-8°C until labeled expiration date.

9) POSITIVE CONTROL SERUM

One vial containing 0.2 ml of anti RNP antibody positive human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

10) NEGATIVE CONTROL SERUM

One vial containing 0.2 ml of anti RNP antibody negative human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

*Those reagents can also be used for REAADS anti-Sm ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, and REAADS anti-CENP-B ELISA Test Kit.

PRECAUTIONS

- (1) This product is for in vitro diagnostic use only.
- (2) Do not use kit components beyond the stated expiration dates.
- (3) Avoid contact of reagents with eyes, skin and clothing. Reagents on skin must be washed away with plenty of water. TMB contains irritant and Stop Solution consists of a 1N sulfuric acid, which is a poison and corrosive.
- (4) Calibrator 2, Positive control and Negative control are derived from human serum, in which HBs antigen, HCV antibody, HIV-1 and HIV-2 antibodies has not been detected. No test method, however, can guarantee the absence of these or any other infectious agents. These reagents and all patient samples should be handled as if they are capable of transmitting AIDS, hepatitis or any other infectious diseases.

- (5) Calibrator 1, Calibrator 2, Positive control, Negative control and Assay Diluent contain sodium azide (0.1%) as a preservative and must be handled with caution - do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (6) Some kit components contain animal origin materials, which are from non-infectious animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (7) Matching lot numbers of Microwell strips, Conjugate and Calibrator 2 must be used together in the assay. Do not substitute reagents from other kits.
- (8) All reagents must be brought to room temperature (20 -25oC) before starting the assay.
- (9) Do not expose the kit to direct sun during assay and storage.
- (10) Avoid microbial and cross contamination of reagents or samples.
- (11) Incubation temperatures above or below normal room temperature (20-25oC), shorter or longer time periods of incubation and inaccurate dilution may give erroneous results.
- (12) The wells must be rinsed with Wash Solution properly enough to avoid false positive.
- (13) Carefully pipette not to foam each sample and reagent to avoid cross contamination between microwells.
- (14) All microwell strips, which are not immediately required, should be returned to the zip lock pouch, which must be carefully resealed to avoid moisture absorption.
- (15) Wash concentrate may become turbid at 2-8oC, which does not cause inconsistent results.
- (16) The RNP antibodies value obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedure.

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader (wavelength: 450 nm, 620 nm/reference)
- Multichannel micropipette (e.g. 100 µl - 300 µl) for dispensing conjugate, substrate, and stop solution.
- Single channel pipette (10 µl & 100 µl)
- Reagent reservoir
- Autowasher or wash bottle
- Deionized or distilled water
- One liter graduated cylinder for preparation of wash solution
- Test tubes for patient sample dilutions (e.g. 1000 µl)
- Disposable pipette tips
- Paper towels
- Basin and disinfectant
- Microplate cover

PROCEDURE

■ PREPARATION OF REAGENTS

- Bring all assay materials to room temperature (20-25°C) prior to use.
- Microwell Strips: Remove required microwell strip from pouch and place them in the frame. Promptly return unused strips to refrigerated storage.
- Wash Solution: Prepare 1:10 dilution of The Wash Concentrate prior to use.(ex. add 100 ml of Wash Concentrate to 900 ml of distilled water). The diluted wash solution is stable for 2 weeks at 2-8°C.
- Do not dilute Calibrator 1, Calibrator 2, Conjugated Reagent, Assay Diluent, Substrate and Stop Solution, which are ready-to-use.

■ PREPARATION OF SAMPLES

- Use fresh patient sera. If storage is needed, they should be frozen below -20°C for up to one month, below -70°C for longer storage. Do not repeat freezing and thawing.
*In case stored below -20°C for more than 6 months or freezing and thawing repeatedly, nonspecific results are obtained because of IgG denaturation.
- Dilute each patient serum 1:101 by adding 10 μl of serum to 1ml of Assay Diluent.
*Diluted samples must be used within a day.
*Diluted samples can be also used for REAADS anti-Sm ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, and REAADS anti-CENP-B ELISA Test Kit.
* Assay Diluent may form precipitate, which does not cause inconsistent results.

■ ASSAY PROCEDURE

STEP 1. (SAMPLE INCUBATION)

Using the multi-channel pipetor, transfer 100 μl of each diluted sample, Positive and Negative Controls into the appropriate microwells of the antigen test plate. Add Calibrators directly to appropriate wells. (Do not dilute Calibrators.)

- * Incubation starts on pipetting to the antigen-coated microwells. Pipetting should be completed as quickly as possible.

Cover wells with a plate sealer and incubate for 60 minutes at room temperature (20-25 $^{\circ}\text{C}$).

STEP 2. (WASHING)

Aspirate or discard the well contents. Fill the well with Wash Solution and then completely aspirate or discard the contents. Wash 4 times. Tap the plate on a paper towel to remove any remaining Wash Solution. When autowasher is used, wash 4 times.

- * Each laboratory is recommended to confirm its own appropriate washing times and set-up.
- * Wash Solution should be used at 20-25 $^{\circ}\text{C}$.

STEP 3. (CONJUGATE INCUBATION)

Pour Conjugated Reagent into a reservoir. Add 100 μl of the Conjugated Reagent to each well with multichannel pipette. Cover wells with the plate sealer and incubate for 60 minutes at room temperature (20-25 $^{\circ}\text{C}$).

STEP 4. (WASHING)

Wash the microplate following the STEP 2 procedure.

STEP 5. (SUBSTRATE INCUBATION)

Pour Substrate into a reservoir and pipette 100 μl of the Substrate to each well with multichannel pipette.

- * The reservoir should be different from the one, which was used for pouring conjugate solution. A new disposable reservoir should be used because Substrate is easily oxidized by metal ion.
- * The Substrate, once poured in a reservoir, should not be returned to the bottle.

Cover wells with the plate sealer and incubate for 30 minutes at room temperature (20-25 $^{\circ}\text{C}$).

STEP 6. (STOP REACTION)

Pour Stop Solution into a reservoir. Pipette 100 μl of the solution to each well with multichannel pipette.

■ READING

Read the absorbance of each well at 450 nm. If a dual wave-length platereader is available, set the test wavelength at 450 nm and the reference at 620 nm.

- * Reading should be done as quickly as possible after stopping the reaction.

- Ensure that the bottom of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading the plate.

■ **CALCULATION OF RESULT**

$$\text{Unit value (U/ml)} = \frac{(A_{450}\langle\text{Sample}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)}{(A_{450}\langle\text{Calibrator 2}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)} \times 100$$

*A₄₅₀ is abbreviation of absorbance value at 450 nm.

*An international reference material for anti RNP antibodies is not available; the assay is calibrated in relative arbitrary units.

■ **QUALITY CONTROL**

Each assay result should meet the following criteria.

A₄₅₀ of Calibrator1: ≤0.100

A₄₅₀ of Calibrator 2: ≥ 0.500

The Positive and Negative Controls must give the following results:

	Positive Control	Negative Control
Anti-RNP value (U/ml)	>60	<15

If any of these are not met, the results are invalid and the test should be repeated. Before repeating assay, check the following procedure.

- Incubation Temperature
- Incubation Period of Time
- Washing

TEST INTERPRETATION AND EXPECTED VALUE

The following is intended only as a guide for interpretation. Each laboratory is recommended to establish its own criteria for test interpretation based on sample populations typically encountered.

Anti-RNP value (U/mL)	Interpretation
< 15	Negative for anti-RNP Ab
Greater than or equal to 22	Positive for anti-RNP Ab

If unit value within 15-22 is calculated, passage of clinical features and results of any other tests should be carefully observed.

■ **LIMITATIONS**

As with other diagnostic test procedures, the results obtained with the READS anti-RNP ELISA Test Kit serve only as an aid to diagnosis and should not be interpreted as diagnostics in themselves.

■ **CLINICAL SPECIFICITY AND SENSITIVITY**

Disease	Positive sample	Positive rate
MCTD	196/199	98.5%
SLE	50/119	42.0%
SjS	4/25	16.0%
Normal Serum	1/98	1.0%

MCTD: Mixed connective tissue disease

SLE: Systemic lupus erythematosus

SjS: Sjögren's syndrome

■ **PRECISION**

Repeatability was demonstrated by testing 3 samples in sextuple. Reproducibility was determined by testing 3 samples on 5 different days (day-to-day) and by testing 3 samples in 3 lots (lot-to-lot). %CV values for reproducibility and repeatability were below 15% for each sample.

■ ASSAY RANGE

The assay range of this kit is from 5 U/ml to 300 U/ml.

■ INTERFERING SUBSTANCES

Hemoglobin (up to 480 mg/dl), Bilirubin (up to 20.0 mg/dl), chyle (up to 2,780 unit as Formazine) and/or Rheumatoid Factor (up to 520 IU/ml) are not affective on the assay result, but avoid using highly hemolysed samples or highly lipemic samples.

REFERENCES

1. Tsay, G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. *Arthritis & Rheum.* 30: 389-396, (1987).
2. Tojo, T. et al.: Diversity of autoantibodies and their clinical significance in systemic autoimmune diseases. *Medical immunology* 14: 761-766, (1987).
3. Sharp, G.C. et al.: Mixed connective tissue disease: an apparently distinct rheumatic disease syndrome associated with a specific antibody to an extractable nuclear antigen (ENA). *Amer. J. Med.* 52: 148-152, (1972).
4. Query, C.C. et al.: A common RNA recognition motif identified within a defined U1 RNA binding domain of the 70K U1 snRNP protein. *Cell* 57: 89-101, (1989).
5. Query, C.C. et al.: A human autoimmune protein associated with U1 RNA contains a region of Homology that is cross-reactive with retroviral p30gag antigen. *Cell* 51: 211-220, (1987).

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BIBLIOGRAFIA

1. Tsay,G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. *Arthritis & Rheum.* 30:389-396 (1987).
2. Tojo,T. et al.: Diversity of autoantibodies and their clinical significance in systemic autoimmune diseases. *Medical Immunology* 14:761-766 (1987).
3. Sharp,G.C. et al.: Mixed connective tissue disease: an apparently distinct rheumatic disease syndrome associated with a specific antibody to an extractable nuclear antigen /ENA). *Am.J.Med.* 52:148-152 (1972).
4. Query C.C. et al.: A common RNA recognition motif identified within a defined U1 RNA binding domain of the 70K U1 snRNP protein. *Cell* 57:89-101 (1989).
5. Query C.C. et al.: A human autoimmune protein associated with U1 RNA contains a region of Homology that is cross-reactive with retroviral p30gag antigen. *Cell* 51:211-220 (1987).


**The following symbols are used on the label.
 **Les symboles suivants ont été utilisés sur l'étiquette
 **Die folgenden Symbole finden Sie auf dem Etikett
 **En las etiquetas se usan los siguientes símbolos
 **I seguenti simboli sono usati sulle etichette.

For in vitro diagnostic use
 Dispositif médical de diagnostic in vitro
 Für in-vitro Diagnostik
 Para uso diagnóstico in vitro
 Per uso diagnostico in vitro

Catalogue number
 Référence du catalogue
 Artikelnummer
 Número de catálogo
 Numero di catalogo

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11575 Main Street, Suite 400

Broomfield, Colorado 80020 USA

REAADS[®] is a registered trademark of Corgenix REAADS, Inc.

© 2006, Corgenix REAADS, Inc.

MT Promedt Consulting GmbH

Altenhofstrasse 80

D-66386 St. Ingbert / Germany

11687 01

Effective: 2006-08-30