# **REAADS anti-Sm**

## **ELISA Test Kit**

Cat. No. 10868: 96 wells  
Cat. No. 10868: 96 puits  
Art.-Nr. 10868: 96 Auftragsstellen  
Cat. No. 10868: 96 pocillos  
Cat. n. 10868: 96 pozzetti

Corgenix, Inc.

11575 Main Street, Suite 400, Broomfield, Colorado 80020, USA

Tel: 303-457-4345 Fax : 303-457-4519 [www.corgenixonline.com](http://www.corgenixonline.com)

- CONTENTS -  
- CONTENU -  
- INHALTSVERZEICHNIS -  
- CONTENIDO -  
- CONTENUTI -

English.....	1
Français.....	7
Deutsch.....	13
Español.....	19
Italiano.....	25
Symbols/Symboles/Symbole/Símbolos /Simboli.....	31

## English

### **INTENDED USE**

The REAADS anti-Sm ELISA Test Kit is semi-quantitative enzyme-linked immunosorbent assay (ELISA) for the detection of anti Sm antibodies in human serum. The REAADS anti-Sm ELISA Test Kit is intended for in vitro diagnostic use as an aid in the determination of certain autoimmune diseases.

### **SUMMARY AND EXPLANATION**

Sm antibodies are ENA antibodies named after “Smith”, a patient who suffered from Systemic Lupus Erythematosus (SLE). Sm antibodies are frequently found with RNP antibodies. It is detected in 10-30% of SLE patients. At high levels, Sm antibody is indicative of SLE and is rarely detected in other diseases. It has been adopted as a marker for diagnosis of SLE by the American Rheumatism Association (ARA). Anti-Sm antibodies are observed at a high titer in the active period of SLE and at a low titer in the non-active period. Raynaud’s phenomena and nephropathy are reported more frequently in SLE patients testing positive for anti-Sm antibodies.

The REAADS anti-Sm ELISA Test Kit measures with high sensitivity, anti-Sm antibodies present in the serum by ELISA.

### **PRINCIPLE**

The REAADS anti-Sm ELISA Test Kit measures anti-Sm antibodies present in the serum by ELISA. Calibrators and patient serum are added to microwell coated with Sm antigens, allowing anti-Sm antibodies to react with the immobilized antigen (Sample incubation). After wash to remove any unbound serum proteins, horseradish peroxidase conjugated anti human IgG, IgA and IgM are added and incubated (Conjugate incubation). Following another washing step, the peroxidase substrate is added and incubated for an additional period of time (Substrate incubation). Acid solution is then added to each well to terminate the enzyme reaction and to stabilize the color development. The assay can be quantified by measuring the reaction photometrically and calculating the results.

### **BRIEF ASSAY PROCEDURE**

<Sample incubation> (20-25 °C) 60 min.	Add 100 µl of diluted sample (1:101) to each well of microwell plate ↓ Wash ↓
<Conjugate incubation> (20-25 °C) 60 min.	Add 100 µl of conjugate solution to each well ↓ Wash ↓
<Substrate incubation> (20-25 °C) 30 min	Add 100 µl of substrate to each well ↓ Add 100 µl of stop solution to each well ↓ Read absorbance ↓ Interpretation of result

## **REAGENTS AND STORAGE**

### **1) Sm MICROWELL STRIPS**

96 wells MICROWELL STRIPS (12 x 8 wells) coated with antigen produced from native purified proteins, the strips packaged in a strip holder and sealed in a foil envelope with desiccant, are stable at 2-8°C until the labeled expiration date.

### **2) CALIBRATOR 1 (0U/ml)\***

Two vials containing 1.5 ml of Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

### **3) CALIBRATOR 2 (100U/ml)**

Two vials containing 1.5 ml each of anti Sm antibody positive human serum with Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

### **4) CONJUGATE REAGENT**

One vial containing 15 ml of horseradish peroxidase conjugated goat anti-human IgG, IgA and IgM (heavy chain specific), HEPES, Proclin 150 and BSA. Ready to use. Stable at 2-8°C until labeled expiration date.

### **5) ASSAY DILUENT\***

Two 50 ml bottles containing PBS, Tween-20 and 0.1% sodium azide. Ready to use. Stable at 2-8°C until labeled expiration date.

### **6) WASH CONCENTRATE (10X) \***

One 100 ml bottle containing PBS and Tween-20 as a 10X concentrate. Stable at 2-8°C until labeled expiration date.

### **7) SUBSTRATE SOLUTION\***

One 20 ml bottle containing 3,3',5,5'-tetramethylbenzidine dihydrochloride/hydrogen peroxide (TMB/H<sub>2</sub>O<sub>2</sub>). Ready to use. Stable at 2-8°C until labeled expiration date.

### **8) STOP SOLUTION\***

One 20 ml bottle containing 1N Sulfuric acid. Ready to use. Stable at 2-8°C until labeled expiration date.

### **9) POSITIVE CONTROL SERUM**

One vial containing 0.2 ml of anti Sm antibody positive human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

### **10) NEGATIVE CONTROL SERUM**

One vial containing 0.2 ml of anti Sm antibody negative human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

\*Those reagents can also be used for REAADS anti-RNP ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, and REAADS anti-CENP-B ELISA Test Kit.

## **PRECAUTIONS**

- (1) This product is for in vitro diagnostic use only.
- (2) Do not use kit components beyond the stated expiration dates.
- (3) Avoid contact of reagents with eyes, skin and clothing. Reagents on skin must be washed away with plenty of water. TMB contains irritant and Stop Solution consists of a 1N sulfuric acid, which is a poison and corrosive.
- (4) Calibrator 2, Positive control and Negative control are derived from human serum, in which HBs antigen, HCV antibody HIV-1 and HIV-2 antibodies has not been detected. No test method, however, can guarantee the absence of these or any other infectious agents. These reagents and all patient samples should be handled as if they are capable of transmitting AIDS, hepatitis or any other infectious diseases.

- (5) Calibrator 1, Calibrator 2, Positive control, Negative control and Assay Diluent contain sodium azide (0.1%) as a preservative and must be handled with caution - do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (6) Some kit components contain animal origin materials, which are from non-infectious animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (7) Matching lot numbers of Microwell strips, Conjugate and Calibrator 2 must be used together in the assay. Do not substitute reagents from other kits.
- (8) All reagents must be brought to room temperature (20 -25 °C) before starting the assay.
- (9) Do not expose the kit to direct sun during assay and storage.
- (10) Avoid microbial and cross contamination of reagents or samples.
- (11) Incubation temperatures above or below normal room temperature (20-25 °C), shorter or longer time periods of incubation and inaccurate dilution may give erroneous results.
- (12) The wells must be rinsed with Wash Solution properly enough to avoid false positive.
- (13) Carefully pipette not to foam each sample and reagent to avoid cross contamination between microwells.
- (14) All microwell strips, which are not immediately required, should be returned to the zip lock pouch, which must be carefully resealed to avoid moisture absorption.
- (15) Wash concentrate may become turbid at 2-8°C, which does not cause inconsistent results.
- (16) The Sm antibodies value obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedure.

## **MATERIALS REQUIRED BUT NOT PROVIDED**

- Microplate reader (wavelength: 450 nm, 620 nm/reference)
- Multichannel micropipette (e.g. 100 µl - 300 µl) for dispensing conjugate, substrate, and stop solution.
- Single channel pipette (10 µl & 100 µl)
- Reagent reservoir
- Autowasher or wash bottle
- Deionized or distilled water
- One liter graduated cylinder for preparation of wash solution
- Test tubes for patient sample dilutions (e.g. 1000 µl)
- Disposable pipette tips
- Paper towels
- Basin and disinfectant
- Microplate cover

## **PROCEDURE**

### **■ PREPARATION OF REAGENTS**

- Bring all assay materials to room temperature (20-25 °C) prior to use.
- Microwell Strips: Remove required microwell strip from pouch and place them in the frame. Promptly return unused strips to refrigerated storage.
- Wash Solution: Prepare 1:10 dilution of The Wash Concentrate prior to use. (ex. add 100 ml of Wash Concentrate to 900 ml of distilled water). The diluted wash solution is stable for 2 weeks at 2 -8°C.
- Do not dilute Calibrator 1, Calibrator 2, Conjugated Reagent, Assay Diluent, Substrate and Stop Solution, which are ready-to-use.

## ■ PREPARATION OF SAMPLES

- Use fresh patient sera. If storage is needed, they should be frozen below -20°C for up to one month, below -70°C for longer storage. Do not repeat freezing and thawing.
  - \*In case stored below -20°C for more than 6 months or freezing and thawing repeatedly, nonspecific results are obtained because of IgG denaturation.
- Dilute each patient serum 1:101 by adding 10µl of serum to 1ml of Assay Diluent.
  - \*Diluted samples must be used within a day.
  - \*Diluted samples can be also used for REAADS anti-RNP ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, and REAADS anti-CENP-B ELISA Test Kit.
  - \*Assay Diluent may form precipitate, which does not cause inconsistent results.

## ■ ASSAY PROCEDURE

### STEP 1. (SAMPLE INCUBATION)

Using the multi-channel pipetor, transfer 100µl of each diluted sample, Positive and Negative Controls into the appropriate microwells of the antigen test plate. Add Calibrators directly to appropriate wells. ( Do not dilute Calibrators.)

\*Incubation starts on pipetting to the antigen-coated microwells. Pipetting should be completed as quickly as possible.

Cover wells with a plate sealer and incubate for 60 minutes at room temperature (20-25 °C).

### STEP 2. (WASHING)

Aspirate or discard the well contents. Fill the well with Wash Solution and then completely aspirate or discard the contents. Wash 4 times. Tap the plate on a paper towel to remove any remaining Wash Solution. When autowasher is used, wash 4 times.

\*Each laboratory is recommended to confirm its own appropriate washing times and set- up.

\*Wash Solution should be used at 20-25 °C.

### STEP 3. (CONJUGATE INCUBATION)

Pour Conjugated Reagent into a reservoir. Add 100 µl of the Conjugated Reagent to each well with multichannel pipette. Cover wells with the plate sealer and incubate for 60 minutes at room temperature (20-25°C).

### STEP 4. (WASHING)

Wash the microplate following the STEP 2 procedure.

### STEP 5. (SUBSTRATE INCUBATION)

Pour Substrate into a reservoir and pipette 100 µl of the Substrate to each well with multichannel pipette.

\*The reservoir should be different from the one, which was used for pouring conjugate solution. A new disposable reservoir should be used because Substrate is easily oxidized by metal ion.

\*The Substrate, once poured in a reservoir, should not be returned to the bottle.

Cover wells with the plate sealer and incubate for 30 minutes at room temperature (20-25°C).

### STEP 6. (STOP REACTION)

Pour Stop Solution into a reservoir. Pipette 100 µl of the solution to each well with multichannel pipette.

## ■ READING

Read the absorbance of each well at 450 nm. If a dual wavelength plate reader is available, set the test wavelength at 450 nm and the reference at 620 nm.

\*Reading should be done as quickly as possible after stopping the reaction.

\*Ensure that the bottom of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading the plate.

**■ CALCULATION OF RESULT**

$$\text{Unit value (U/ml)} = \frac{(A_{450}\langle\text{Sample}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)}{(A_{450}\langle\text{Calibrator 2}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)} \times 100$$

\*A<sub>450</sub> is abbreviation of absorbance value at 450 nm.

\*An international reference material for anti Sm antibodies is not available, the assay is calibrated in relative arbitrary units.

**■ QUALITY CONTROL**

Each assay result should meet the following criteria.

A<sub>450</sub> of Calibrator1: ≤0.100

A<sub>450</sub> of Calibrator2: ≥0.500

The Positive and Negative Controls must give the following results:

	Positive Control	Negative Control
Anti-Sm value (U/ml)	> 100	<7

If any of these are not met, the results are invalid and the test should be repeated.

Before repeating assay, check the following procedure.

- Incubation Temperature
- Incubation Period of Time
- Washing

**TEST INTERPRETATION AND EXPECTED VALUE**

The following is intended only as a guide for interpretation. Each laboratory is recommended to establish its own criteria for test interpretation based on sample populations typically encountered.

Anti-Sm value (U/ml)	Interpretation
< 7	Negative for anti-Sm Ab
Greater than or equal to 30	Positive for anti-Sm Ab

If unit value within 7-30 is calculated, passage of clinical features and results of any other tests should be carefully observed.

**■ LIMITATIONS**

As with other diagnostic test procedures, the results obtained with the REAADS anti-Sm ELISA Test Kit serve only as an aid to diagnosis and should not be interpreted as diagnostics in themselves.

**PERFORMANCE CHARACTERISTICS**

**■ CLINICAL SPECIFICITY AND SENSITIVITY**

Disease	Positive sample	Positive rate
MCTD	6/38	15.8%
SLE	63/260	24.2%
SjS	6/130	4.6%
SSc	3/68	4.4%
Normal Serum	0/160	0.0%

MCTD: Mixed connective tissue disease

SLE: Systemic lupus erythematosus

SjS: Sjögren’s syndrome

SSc: Systemic sclerosis

## ■ PRECISION

Repeatability was demonstrated by testing 3 samples in sextuple. Reproducibility was determined by testing 3 samples on 5 different days (day-to-day) and by testing 3 samples in 3 lots (lot-to-lot). %CV values for reproducibility and repeatability were below 15% for each sample.

## ■ ASSAY RANGE

The assay range of this kit is from 5U/ml to 300U/ml.

## ■ INTERFERING SUBSTANCES

Hemoglobin (up to 480 mg/dl), Bilirubin (up to 20.0mg/dl), chyle (up to 2,780 unit as Formazine) and/or Rheumatoid Factor (up to 520 IU/ml) are not affective on the assay result, but avoid using highly hemolysed samples or highly lipemic samples.

## REFERENCES

1. Paxton, H. et al: Evaluation of the Rheumastrip ANA profile test: a rapid test procedure for simultaneously determining antibodies to autoantigens U1-Ribonucleoprotein (U1-RNP), Sm, SS-A/Ro, SS-B/La, and to native DNA. Clin. Chem. 36: 792-797, (1990).
2. Eng M. Tan: Diagnostic markers for autoimmune diseases and probes for cell biology. Advances in immunology 44: 93-151, (1989).
3. Tsay, G. J. et al: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. Arthritis & Rheum. 30: 389-396, (1987).

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## BIBLIOGRAFIA

1. Paxton,H., et al. : Evaluation of the RheumaStrip ANA profile test : a rapid test strip procedure for simultaneously determining antibodies to autoantigens U1-Ribonucleoprotein (U1-RNP), Sm, SS-A/Ro, SS-B/La, and to Native DNA. Clin.Chem. 36/5:792-797,(1990).
2. Eng,M. Tan: Diagnostic markers for autoimmune diseases and probes for cell biology. Advances in Immunology 44:93-151 (1989).
3. Tsay,G.J. et al.: An immunoassay differentiating sera with antibodies to Sm alone, antibodies to Sm/RNP complex, and antibodies to RNP alone. Arthritis & Rheum. 30:389-396 (1987).

\*\*The following symbols are used on the label.  
 \*\*Les symboles suivants ont été utilisés sur l'étiquette.  
 \*\*Die folgenden Symbole finden Sie auf dem Etikett.  
 \*\*En las etiquetas se usan los siguientes símbolos  
 \*\*I seguenti simboli sono usati sulle etichette.

For in vitro diagnostic use  
 Dispositif médical de diagnostic in vitro  
 Für in-vitro Diagnostik  
 Para uso diagnóstico in vitro  
 Per uso diagnostico in vitro

Catalogue number  
 Référence du catalogue  
 Artikelnummer  
 Número de catálogo  
 Numero di catalogo

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© 2002, Corgenix REAADS, Inc.
REAADS® is a registered trademark of Corgenix REAADS, Inc.
**MT Promedt Consulting GmbH**