REAADS anti-PR-3

ELISA Test Kit

Cat. No. 10881: 96 wells
Cat. No. 10881: 96 puits
Art.-Nr. 10881: 96 Auftragsstellen
Cat. No. 10881: 96 pocillos
Cat. n. 10881: 96 pozzetti

11575 Main Street, Suite 400, Broomfield, Colorado 80020 USA
Tel: 303-457-4345 Fax : 303-457-4519 www.corgenixonline.com

- CONTENTS -
- CONTENU -
- INHALTSVERZEICHNIS -
- CONTENIDO -
- CONTENUTI -

English.....	3
Français.....	9
Deutsch.....	15
Español.....	21
Italiano.....	27
Symbols/Symboles/Symbole/Símbolos /Simboli...	33

ENGLISH

INTENDED USE

The REAADS anti-PR-3 ELISA Test Kit is a semi-quantitative enzyme-linked immunosorbent assay (ELISA) for the detection of IgG anti-proteinase (PR-3) antibodies in human serum. The REAADS anti-PR-3 ELISA Test Kit is intended for in vitro diagnostic use as an aid in the diagnosis of certain systemic vasculitides such as Wegener's granulomatosis.

SUMMARY AND EXPLANATION

In 1982 Davis et al. described the presence of antineutrophil cytoplasmic antibodies (ANCA) in serum of patients with necrotising glomerulonephritis. The use of indirect immunofluorescence (IIF) identified two distinct pattern types described in 1988 by Falk et al, pANCA (perinuclear) and cANCA (cytoplasmic), in renal patients with systemic vasculitis. The major pANCA target antigen is myeloperoxidase and the major cANCA target antigen is proteinase III (PR-3). Several other antigens (lactoferrin, elastase, cathepsin G) are associated with ANCA reactivity to lesser degrees. Although the use of an ELISA can identify the presence of antibodies specific for MPO and PR-3, it is still recommended to use an IIF method for patient screening. Testing positive IIF samples by ELISA can provide useful additional information. ELISA assays that identify the presence of anti-MPO and anti-PR-3 antibodies can be used not only to confirm IIF positive patterns, but also distinguish specific systemic vasculitides. Anti-PR-3 antibodies are detected in the majority of active Wegener's granulomatosis patients and only occasionally positive for anti-MPO. Diseases such as microscopic polyarteritis and crescentic glomerulonephritis are more likely associated with the presence of anti-MPO antibodies.

PRINCIPLE

The REAADS anti-PR-3 ELISA Test Kit measures anti-PR-3 antibodies present in the serum by ELISA. Calibrators and patient serum are added to microwells coated with PR-3 antigens, allowing anti-PR-3 antibodies to react with the immobilized antigen (Sample incubation). After wash to remove any unbound serum proteins, horseradish peroxidase conjugated anti-human IgG are added and incubated (Conjugate incubation). Following another washing step, the peroxidase substrate is added and incubated for an additional period of time (Substrate incubation). Acid solution is then added to each well to terminate the enzyme reaction and to stabilize the color development. The assay can be quantified by measuring the reaction photometrically and plotting the results.

BRIEF ASSAY PROCEDURE

	Add 100 µl of diluted sample (1:101) to each well of microwell plate
	↓
<Sample incubation> (20-25 °C) 60 min.	Wash
	↓
	Add 100 µl of conjugate solution to each well
	↓
<Conjugate incubation> (20-25 °C) 60 min.	Wash
	↓
	Add 100 µl of substrate to each well
	↓
<Substrate incubation> (20-25 °C) 30 min	Add 100 µl of stop solution to each well
	↓
	Read absorbance
	↓
	Interpretation of result

REAGENTS AND STORAGE

1) PR-3 MICROWELL STRIPS

96 wells MICROWELL STRIPS (12 x 8 wells) coated with antigen produced from native purified proteins, the strips packaged in a strip holder and sealed in a foil envelope with desiccant, are stable at 2-8°C until the labeled expiration date.

2) CALIBRATOR 1 (0 U/ml)*

Two vials containing 1.5 ml of Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

3) CALIBRATOR 2 (100 U/ml)

Two vials containing 1.5 ml each of anti-PR-3 antibody positive human serum with Assay Diluent including 0.1% sodium azide. Ready to use, make no further dilution. Stable at 2-8°C until labeled expiration date.

4) CONJUGATE REAGENT

One vial containing 15 ml of horseradish peroxidase conjugated goat anti-human IgG (heavy chain specific), HEPES, Proclin 150 and BSA. Ready to use. Stable at 2-8°C until labeled expiration date.

5) ASSAY DILUENT*

Two 50 ml bottles containing PBS, Tween-20 and 0.1% sodium azide. Ready to use. Stable at 2-8°C until labeled expiration date.

6) WASH CONCENTRATE (10X)*

One 100 ml bottle containing PBS and Tween-20 as a 10X concentrate. Stable at 2-8°C until labeled expiration date.

7) SUBSTRATE SOLUTION*

One 20 ml bottle containing 3,3',5,5'-tetramethylbenzidine dihydrochloride/hydrogen peroxide (TMB/H₂O₂). Ready to use. Stable at 2-8°C until labeled expiration date.

8) STOP SOLUTION*

One 20 ml bottle containing 1N Sulfuric acid. Ready to use. Stable at 2-8°C until labeled expiration date.

9) POSITIVE CONTROL SERUM

One vial containing 0.2 ml of anti-PR-3 antibody positive human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

10) NEGATIVE CONTROL SERUM

One vial containing 0.2 ml of anti-PR-3 antibody negative human serum with 0.1% sodium azide. Stable at 2-8°C until labeled expiration date.

*Those reagents can also be used for REAADS anti-RNP ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, REAADS anti-CENP-B ELISA Test Kit, and REAADS anti-MPO ELISA Test Kit.

PRECAUTIONS

- (1) This product is for in vitro diagnostic use only.
- (2) Do not use kit components beyond the stated expiration dates.
- (3) Avoid contact of reagents with eyes, skin and clothing. Reagents on skin must be washed away with plenty of water. TMB contains irritant and Stop Solution consists of a 1N sulfuric acid.
- (4) Calibrator-2, Positive control and Negative control are derived from human serum, in which HBs antigen, HCV antibody HIV-1 and HIV-2 antibodies have not been detected. No test method, however, can guarantee the absence of these or any other infectious agents. These reagents and all patient samples should be handled as if they are capable of transmitting AIDS, hepatitis or any other infectious diseases.

- (5) Calibrator 1, Calibrator 2, Positive control, Negative control and Assay Diluent contain sodium azide (0.1%) as a preservative and must be handled with caution - do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (6) Some kit components contain animal origin materials, which are from non-infectious animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (7) Matching lot numbers of Microwell strips, Conjugate and Calibrator 2 must be used together in the assay. Do not substitute reagents from other kits.
- (8) All reagents must be brought to room temperature (20-25°C) before starting the assay.
- (9) Do not expose the kit to direct sun during assay and storage.
- (10) Avoid microbial and cross contamination of reagents or samples.
- (11) Incubation temperatures above or below normal room temperature (20-25°C), shorter or longer time periods of incubation and inaccurate dilution may give erroneous results.
- (12) The wells must be rinsed with Wash Solution properly enough to avoid false positive.
- (13) Carefully pipette not to foam each sample and reagent to avoid cross contamination between microwells.
- (14) All microwell strips, which are not immediately required, should be returned to the zip lock pouch, which must be carefully resealed to avoid moisture absorption.
- (15) Wash concentrate may become turbid at 2-8°C, which does not cause inconsistent results.
- (16) The PR-3 antibody values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedure.

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader (wavelength: 450 nm, 620 nm/reference)
- Multichannel micropipette (e.g. 100 µl - 300 µl) for dispensing conjugate, substrate, and stop solution.
- Single channel pipette (10 µl & 100 µl)
- Reagent reservoir
- Autowasher or wash bottle
- Deionized or distilled water
- One liter graduated cylinder for preparation of wash solution
- Test tubes for patient sample dilutions (e.g. 1000 µl)
- Disposable pipette tips
- Paper towels
- Basin and disinfectant
- Microplate cover

PROCEDURE

■ PREPARATION OF REAGENTS

- Bring all assay materials to room temperature (20-25°C) prior to use.
- Microwell Strips: Remove required microwell strip from pouch and place them in the frame. Promptly return unused strips to refrigerated storage.
- Wash Solution: Prepare 1:10 dilution of the Wash Concentrate prior to use. (ex. add 100 ml of Wash Concentrate to 900 ml of distilled water). The diluted wash solution is stable for 2 weeks at 2-8°C.
- Do not dilute Calibrator-1, Calibrator-2, Conjugated Reagent, Assay Diluent, Substrate and Stop Solution, which are ready-to-use.

■ PREPARATION OF SAMPLES

- Use fresh patient sera. If storage is needed, they should be frozen below -20°C for up to one month, below -70°C for longer storage. Do not repeat freezing and thawing.
 - *In case stored below -20°C for more than 6 months or freezing and thawing repeatedly, nonspecific results are obtained because of IgG denaturation.
- Dilute each patient serum 1:101 by adding 10µl of serum to 1ml of Assay Diluent.
 - *Diluted samples must be used within a day.
 - *Diluted samples can also be used for REAADS anti-RNP ELISA Test Kit, REAADS anti-SS-A ELISA Test Kit, REAADS anti-SS-B ELISA Test Kit, REAADS anti-Scl-70 ELISA Test Kit, REAADS anti-Jo-1 ELISA Test Kit, REAADS anti-CENP-B ELISA Test Kit, and REAADS anti-MPO ELISA Test Kit.
 - *Assay Diluent may form precipitate, which does not cause inconsistent results.

■ ASSAY PROCEDURE

STEP 1. (SAMPLE INCUBATION)

Using the multi-channel pipetor, transfer 100µl of each diluted sample, Positive and Negative Controls into the appropriate microwells of the antigen test plate. Add Calibrators directly to appropriate wells. (Do not dilute Calibrators.)

*Incubation starts on pipetting to the antigen-coated microwells. Pipetting should be completed as quickly as possible.

Cover wells with a plate sealer and incubate for 60 minutes at room temperature (20-25°C).

STEP 2. (WASHING)

Aspirate or discard the well contents. Fill the well with Wash Solution and then completely aspirate or discard the contents. Wash 4 times. Tap the plate on a paper towel to remove any remaining Wash Solution. When autowasher is used, wash 4 times.

* Each laboratory is recommended to confirm its own appropriate washing times and set- up.

Wash Solution should be used at 20-25°C.

STEP 3. (CONJUGATE INCUBATION)

Pour Conjugated Reagent into a reservoir. Add 100 µl of the Conjugated Reagent to each well with multichannel pipette. Cover wells with the plate sealer and incubate for 60 minutes at room temperature (20-25°C).

STEP 4. (WASHING)

Wash the microplate following the STEP 2 procedure.

STEP 5. (SUBSTRATE INCUBATION)

Pour Substrate into a reservoir and pipette 100 µl of the Substrate to each well with multichannel pipette.

* The reservoir should be different from the one, which was used for pouring conjugate solution. A new disposable reservoir should be used because Substrate is easily oxidized by metal ions.

* The Substrate, once poured in a reservoir, should not be returned to the bottle.

Cover wells with the plate sealer and incubate for 30 minutes at room temperature (20-25°C).

STEP 6. (STOP REACTION)

Pour Stop Solution into a reservoir. Pipette 100 µl of the solution to each well with multichannel pipette.

■ READING

Read the absorbance of each well at 450 nm. If a dual wavelength plate reader is available, set the test wavelength at 450 nm and the reference at 620 nm.

*Reading should be done as quickly as possible after stopping the reaction.

*Ensure that the bottom of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading the plate.

■ **CALCULATION OF RESULT**

$$\text{Unit value (U/ml)} = \frac{(A_{450}\langle\text{Sample}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)}{(A_{450}\langle\text{Calibrator 2}\rangle - A_{450}\langle\text{Calibrator 1}\rangle)} \times 100$$

*A₄₅₀ is abbreviation of absorbance value at 450 nm.

*An international reference material for anti PR-3 antibodies is not available. The assay is calibrated in relative arbitrary units.

■ **QUALITY CONTROL**

Each assay result should meet the following criteria.

A₄₅₀ of Calibrator1: ≤0.100

A₄₅₀ of Calibrator2: ≥0.500

The Positive and Negative Controls must give the following results:

	Positive Control	Negative Control
Anti-PR-3 value (U/ml)	> 50	< 23

If any of these are not met, the results are invalid and the test should be repeated. Before repeating assay, check the following procedure.

- Incubation Temperature
- Incubation Period of Time
- Washing

TEST INTERPRETATION AND EXPECTED VALUE

The following is intended only as a guide for interpretation. Each laboratory is recommended to establish its own criteria for test interpretation based on sample populations typically encountered.

Anti-PR-3 value (U/mL)	Interpretation
<23	Negative for anti-PR-3 Ab
Greater than or equal to 23	Positive for anti-PR-3 Ab

■ **LIMITATIONS**

As with other diagnostic test procedures, the results obtained with the READS anti-PR-3 ELISA Test Kit serve only as an aid to diagnosis and should not be interpreted as diagnostics in themselves.

PERFORMANCE CHARACTERISTICS

■ **CLINICAL SPECIFICITY AND SENSITIVITY**

Disease	Positive sample	Positive rate
Vasculitis cANCA+	27/31	87%
Vasculitis cANCA-	0/47	0%
MCTD	0/10	0%
SLE	1/10	10%
SjS	0/10	0%
PM/DM	0/10	0%
RA	0/9	0%
Normal Serum	0/80	0%

MCTD: Mixed connective tissue disease

SLE: Systemic lupus erythematosus

SjS: Sjögren's syndrome

PM/DM: Polymyositis/Dermatomyositis

RA: Rheumatoid Arthritis

■ PRECISION

Repeatability was demonstrated by testing 3 samples in replicate. Reproducibility was determined by testing 3 samples on 7 different days (day-to-day) and by testing 3 samples in 3 lots (lot-to-lot). %CV values for reproducibility and repeatability were below 15% for each sample.

■ ASSAY RANGE

The assay range of this kit is from 5 U/ml to 200 U/ml.

■ INTERFERING SUBSTANCES

Hemoglobin (up to 440 mg/dl), Bilirubin (up to 19.5 mg/dl), chyle (up to 2,350 unit as Formazine) and/or Rheumatoid Factor (up to 500 IU/ml) are not affective on the assay result, but avoid using highly hemolysed samples or highly lipemic samples.

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© 2007, Corgenix, Inc.
REAADS® is a registered trademark of Corgenix, Inc.
MT Promedt Consulting GmbH