REAADS® Anti-dsDNA Quantitative Test Kit
For In Vitro Diagnostic Use

An enzyme-linked immunosorbent assay (ELISA) for the quantitative determination of IgG/IgM anti-dsDNA antibodies in human serum.

INTENDED USE

Detection and quantitation of IgG/IgM anti-double stranded deoxyribonucleic acid (dsDNA) antibodies as an aid in the diagnosis and management of systemic lupus erythematosus (SLE).

SUMMARY AND EXPLANATION OF THE ANTI-dsDNA TEST

Antibodies directed against deoxyribonucleic acid (DNA) were first detected in the serum of patients with systemic lupus erythematosus (SLE) in the late 1950s.^{1,2,3,4} The presence of anti-dsDNA autoantibodies is one of the four highly specific serologic markers included in the 1982 ARA (ACR) revised criteria for the classification of SLE.⁵ The detection of elevated levels of anti-dsDNA antibodies and decreased serum levels of complement component C3 have been found to be 100% specific for SLE.⁶

Antibodies against native, or double stranded DNA (dsDNA) are found almost exclusively in patients with SLE, while antibodies against denatured, or single stranded DNA (ssDNA) can also be detected in patients with SLE as well as in patients with many other diseases, including rheumatic diseases and chronic infections.^{5,7,8}

The serum level of anti-dsDNA antibodies in patients with SLE correlates significantly with the level of disease activity, particularly when there is renal involvement.^{9,10,11} Therefore, the test for anti-dsDNA antibodies is useful for monitoring disease activity and the progress of therapy in patients with SLE.

A variety of commercially available tests incorporate different technologies for the detection of anti-dsDNA antibodies, including solution phase assays that use radiolabeled DNA, and indirect immunofluorescent assays in which anti-dsDNA antibodies bind to the DNA of the hemoflagellate *Crithidia luciliae*. Most of these techniques have inherent limitations. The radiolabeling process of the solution phase assays can denature segments of the DNA, reducing the specificity of these assays for anti-dsDNA antibodies.¹¹ The immuno-fluorescent *Crithidia* staining assay is labor intensive, and requires the subjective interpretation of the staining pattern.

PRINCIPLE OF THE TEST

The test is performed as an indirect ELISA. Diluted serum samples, calibrators, and controls are incubated in dsDNA coated microwells, allowing anti-dsDNA antibodies present in the samples to react with the immobilized antigen. After the removal of unbound serum proteins by washing, antibodies specific for human IgG and IgM labeled with horseradish peroxidase (HRP) are added forming complexes with dsDNA bound antibodies. Following another washing step, the bound enzyme-antibody conjugate is assayed by the addition of a single solution containing tetramethylbenzidine (TMB) and hydrogen peroxide (H₂O₂) as the chromogenic substrate. Color develops in the wells at an intensity proportional to the serum concentration of anti-dsDNA antibodies.

Results are obtained by reading the O.D. (optical density or absorbance) of each well in a spectrophotometer. A calibrator serum is provided, with the anti-dsDNA antibody concentration expressed in either AU/mL (arbitrary units that have been standardized against a Centers for Disease Control reference preparation) or in IU/mL (international units that have been standardized against the World Health Organization reference preparation). Dividing the concentration value of the calibrator by the absorbance value of that calibrator provides a conversion factor. The O.D. values of the controls and patient samples are multiplied by the conversion factor to obtain antibody concentration values, expressed in either AU/MI or IU/mL.

REAGENTS

Store at 2-8°C. Do Not Freeze.

Each REAADS Anti-dsDNA 96-microwell Test Kit contains the following reagents (**volumes may vary depending on kit size and configuration**):

- 12 1X8 stabilized dsDNA (from calf thymus) coated microwell strips with frame
- 60 mL Sample Diluent I (green solution); contains bovine calf serum*
- 0.250 mL Calibrator Serum; see vial label for anti-dsDNA antibody concentration in AU/mL and IU/mL*
- 0.250 mL Positive Control Serum; see vial label for anti-dsDNA antibody concentration range in AU/mL and IU/mL*
- 0.250 mL Normal Control Serum; see vial label for anti-dsDNA antibody concentration range in AU/mL and IU/mL*
- 15 mL anti-human IgG/IgM (goat) HRP-Conjugated Antibody Solution (red solution); contains 0.01% thimerosal and gentamycin sulfate as preservatives
- 15 mL One Component Substrate (TMB and H₂O₂); ready to use
- 15 mL Stopping Solution (0.36 N sulfuric acid)
- 2 x 30 mL Wash Concentrate (33X PBS)

*** CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use

1. Human source material used to prepare the calibrator and controls included in this kit has been tested and shown to be negative for antibodies to HBsAg, HCV and HIV I & 2 by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Certain components of this product contain sodium azide as a preservative. Sodium azide has been reported to form lead and copper azides when left in contact with these metals. These metal azides are explosive. Any solutions containing azide must be thoroughly flushed with copious amounts of water to prevent the build-up of explosive metal azides in the plumbing system.
6. One component substrate can cause irritation to the eyes and skin. Absorption through the skin is possible. Use gloves when handling substrate and wash thoroughly after handling. Keep reagent away from ignition sources. Avoid contact with oxidizing agents.
7. Certain components are labeled with the following:
Irritating to eyes (R 36). Avoid contact with skin (S 24). Avoid contact with eyes (S 25). In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S 26). If swallowed, seek medical advice immediately and show this container or label (S 46).

Irritant . Biological Risk .

SPECIMEN COLLECTION AND PREPARATION

Serum is the preferred sample matrix. Blood should be collected by venipuncture, and the serum separated from the cells by centrifugation after clot formation. If not tested immediately, specimens should be stored at 2 to 8°C. If specimens are to be stored for more than 72 hours, they should be frozen at -20°C or below. Avoid repeated freezing and thawing. Do not use hemolyzed, icteric, or lipemic serum or plasma as these conditions may cause aberrant results. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

INSTRUCTIONS FOR USE

Materials Provided:

REAADS Anti-dsDNA Test Kit; see “Reagents” for a complete listing.

Materials Required but not Provided:

- Reagent grade water to prepare PBS wash solution (1L) and to zero or blank the plate reader during the final assay step
- Graduated cylinders
- Precision pipettors capable of delivering between 5 μ L and 1000 μ L, with appropriate tips
- Miscellaneous glassware appropriate for small volume handling
- Flask or bottle, 1 liter
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semi-automated washing system
- Disposable gloves
- Plate reading spectrophotometer capable of reading absorbance at 450 nm (with a 650 nm reference if available)
- Multichannel pipettors capable of delivering to 8 wells simultaneously

Procedural Notes

1. Bring serum samples and kit reagents to room temperature (18-26°C) and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage as soon as possible.
2. All dilutions of calibrator, control, and test sera must be made just prior to use in the assay.
3. A single water blank well can be set up on each plate with each run. No sample or kit reagents are to be added to this well. Instead, add 200 μ L of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to “zero” or “blank” against air or a water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution from a plastic squeeze bottle with a wide tip into the bottom of the microwells. No interference is caused should a water blank well receive PBS. An automated microtiter plate washing system can also be used.
5. **IMPORTANT:** Failure to adequately remove residual PBS can cause inconsistent color development of the substrate solution.
6. Use a multichannel pipettor capable of delivering to 8 wells simultaneously when possible. This speeds the process and provides more uniform incubation and reaction times for all wells.
7. Carefully controlled timing of all steps is critical. All calibrators, controls, and samples must be added within a five minute period. Batch size of samples should not be larger than the amount that can be added within this time period.
8. For all incubations, the start of the incubation period begins with the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below normal room temperature (18-26°C) may contribute to inaccurate results.
11. Avoid microbial and cross-contamination of reagents when opening and removing aliquots from the primary vials.
12. Do not use Tween 20 or other detergents in this assay.
13. Do not use kit components beyond expiration date.
14. Do not use kit components from different kit lot numbers.

Reagent Preparation

Wash Solution (PBS): Measure 30 mL of Wash Concentrate (33X PBS) and dilute to 1 liter with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused PBS solution in the refrigerator at 2-8°C. Discard if the solution shows signs of microbial or cross contamination.

Assay Procedure

1. One well should be run as a reagent blank control; Sample Diluent I without serum is added to the well. This well will be treated the same as sample wells in subsequent assay steps. The remaining wells are sufficient to test 46 patient samples run in duplicate wells or 92 patient samples run in single wells.
2. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
3. Prepare a 1:50 dilution of the calibrator, controls, and patient samples in Sample Diluent I (green solution); e.g. 10 μ L sample added to 490 μ L Sample Diluent I equals a 1:50 sample dilution.
4. Add 100 μ L of each diluted calibrator, control, and patient sample to the appropriate microwells.
5. Add 100 μ L of Sample Diluent I to the reagent blank well.
6. Incubate 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and dump the sample fluid. Do not allow samples to contaminate other microwells.
7. Wash 4 times with PBS. Each well should be filled with PBS per wash. Invert microwells between each wash to empty fluid. Use a snapping motion of the wrist to shake the liquid from the wells. The frame must be squeezed at the center on the top and bottom to retain microwell modules during washing. Blot on absorbent paper to remove residual wash fluid after the final wash. Do not allow the wells to dry out between steps.
8. Add 100 μ L IgG/IgM HRP-Conjugated to the calibrator, controls, reagent blank and patient sample wells.
9. Incubate for 15 minutes at room temperature. After the incubation is complete, carefully invert microwells to dump the enzyme-antibody conjugate solution.
10. Wash 4 times with PBS as in step 7. Use a snapping motion to drain the liquid, and blot on absorbent paper after the final wash. Do not allow the wells to dry out.
11. Add 100 μ L Substrate to each well and incubate for 10 minutes at room temperature. Add the substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
12. Add 100 μ L of the Stopping Solution (0.36 N sulfuric acid) to each well to stop the enzyme reaction. Be sure to add the Stopping Solution to the wells in the same order and at the same rate as the Substrate was added. Blue Substrate will turn yellow and colorless solution will remain colorless. Blank or zero the plate reader against air or a water blank well. Read the O.D. of each well at 450 nm (and 650 nm reference if dual beam). The O.D. values should be measured within 30 minutes after the addition of Stopping Solution.

Results

1. Calculate the mean O.D. values if duplicates of the calibrator, controls, and patient samples were performed.
2. Divide the concentration value (either AU/mL or IU/mL) of the calibrator serum (printed on the vial label) by the O.D. or mean O.D. of the calibrator serum to obtain the conversion factor.
3. Multiply the O.D. or mean O.D. for each of the controls and patient samples by the conversion factor to obtain a concentration value in AU/mL or IU/mL.

$$\begin{aligned} \text{Conversion Factor} &= \frac{\text{Anti-dsDNA concentration of the calibrator}}{\text{Absorbance value of the calibrator (O.D.)}} \\ \text{Anti-dsDNA concentration of sample} &= \text{Conversion Factor} \times \text{Absorbance of the sample (O.D.)} \end{aligned}$$

- The conversion factor must be calculated for the calibrator for each assay run. Using a conversion factor from another assay will invalidate the results.
- Samples with anti-dsDNA values greater than 100 AU/mL may be reported as “greater than 100 AU/mL” or greater than 450 IU/mL may be reported as “greater than 450 IU/mL”.

Quality Control

- The O.D. value of the calibrator should be at least 0.400 to assure that the kit is functioning properly. Calibrator O.D. readings of less than 0.400 may indicate that the kit is no longer suitable for use.
- The O.D. of the reagent blank should be ≤ 0.100 when the spectrophotometer has been blanked against the water well. Readings greater than 0.100 may indicate possible reagent contamination or inadequate plate washing.
- The anti-dsDNA values obtained for the control sera should be within the ranges indicated on the container labels. Occasional small deviations outside these ranges are acceptable.
- O.D. values for duplicates, if run, of the controls or patient samples should be within 20% of the mean O.D. value for samples with absorbance readings greater than 0.200.
- Each laboratory should periodically determine their own normal cut-off values for the appropriate population of patients. See Expected Values, as an example.
- Assure that all quality control parameters have been met before reporting test results.

EXPECTED VALUES

Normal Range: Less than 26 AU/mL
or
Less than 117 IU/mL

The normal cut-off was defined as the mean AU/mL or IU/mL value plus two standard deviations from testing a population of healthy blood donors serum samples.

PERFORMANCE CHARACTERISTICS

NOTE: All performance characteristics were determined by using duplicate well determinations.

Clinical Specificity

Serum samples from 90 healthy blood donors were assayed for the presence of IgG and IgM anti-dsDNA antibodies. The clinical specificity of the test was 98%.

Clinical Sensitivity

Serum samples from 100 individuals diagnosed as having SLE were assayed. Patients with both active and inactive disease were included. Sixty-five (65%) of these samples were positive for anti-dsDNA antibodies. The expected prevalence of elevated levels of anti-dsDNA antibodies in this group (SLE with varying levels of disease activity) is 50-75%.⁴

Precision

Two (2) samples with known anti-dsDNA antibody concentration values (one low and one high) were assayed in 32 replicates on three different occasions. The intraassay and interassay coefficients of variation (CV) are presented in the table below. The reported intraassay CV is the mean of the three separate intraassay CVs.

<u>Sample</u>	<u>Mean Intraassay CV</u>	<u>Interassay CV</u>
Low (25.8 AU/mL or 116.1 IU/mL)	8.6%	8.8%
High (52.1 AU/mL or 234.5 IU/mL)	6.4%	9.8%

LIMITATIONS OF THE TEST

The anti-dsDNA concentration values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures. An elevated serum level of anti-dsDNA antibodies is only one of 11 criteria to be considered for a diagnosis of SLE.¹ Some apparently normal individuals and individuals with other rheumatic diseases may produce measurable amounts of anti-dsDNA antibodies; however, the levels are generally much lower than in patients with SLE.⁸

Warranty

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone +1-303-457-4345, fax +1-303-457-4519 or contact a Corgenix authorized distributor.

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REFERENCES

1. Miescher P, Strässle R. New serological methods for the detection of LE factor. *Vox Sang* 2:283, 1957.
2. Ceppellini R, Polli E, Celada F. A DNA-reacting factor in serum of a patient with lupus erythematosus diffusus. *Proc Soc Exp Biol Med* 96:572, 1957.
3. Robbins WC, Holman HR, Deicher H, Kunkel HG. Complement fixation with cell nuclei and DNA in lupus erythematosus. *Proc Soc Exp Biol Med* 96:575, 1957.
4. Seligmann M. Mise en évidence dans le sérum de malades atteints de lupus érythémateux disséminé d'une substance déterminant une réaction de précipitation avec l'acide désoxyribonucléique. *CR Acad Sci (Paris)* 245:243, 1957.
5. Tan EM. The 1982 revised criteria for the classification of systemic lupus erythematosus. *Arthritis Rheum* 25(11):1271, 1982.
6. Weinstein A, Bordwell B, Stone B, Tibbetts C, Rothfield NF. Antibodies to native DNA and serum complement (C3) levels: application to diagnosis and classification of systemic lupus erythematosus. *Am J Med* 74:206, 1983.
7. Glass D, Shur P. Autoimmunity and Systemic Lupus Erythematosus. *Autoimmunity: Genetic, Immunologic, Virologic and Clinical Aspects. N Talal, M.D., ed* (Academic Press: New York) 531, 1977.
8. Schoenfeld Y, André-Schwartz J, Stollar BD, Schwartz RS. Anti-DNA Antibodies. *Systemic Lupus Erythematosus. RG Lahita, M.D. PhD., ed* (John Wiley and Sons: New York) 213, 1987.
9. Tzioufas AG, Manoussakis MN, Drosos AA, Silis G, Gharavi AE, Moutsopoulos HM. Enzyme immunoassays for the detection of IgG and IgM anti-dsDNA antibodies: clinical significance and specificity. *Clin Exp Rheum* 5:247, 1987.
10. Swaak AJG, Groenwold J, Bronsfield W. Predictive value of complement profiles and anti-dsDNA in systemic lupus erythematosus. *Ann Rheum Dis* 45:359, 1986.
11. Wilson M. Antinuclear Antibodies and Anticytoplasmic Antibodies in Lupus Erythematosus. *Dubois' Lupus Erythematosus, 3rd edition DJ Wallace and EL Dubois, ed* (Lea and Febiger: Philadelphia) 227, 1987.
12. Tan EM, Fritzler MJ, Hardin JA, McDougal JS, et al. Comprehensive Report: Reference Sera for Antinuclear Antibodies I. Antibodies to Native DNA, Sm, Nuclear RNP, SS-B/La, SS-A/Ro, Scl-70, Nucleoli, and Centromere/Kinetochore. Committee on Antinuclear Antibody Serology, The Arthritis Foundation, Atlanta, GA, 1984.

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Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020, USA

READS® is a registered trademark of Corgenix, Inc.
 © 2005, Corgenix, Inc.
 US Patent # 5,183,735

13022901 14
 Effective: 2006-10-30

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MT Promedt Consulting GmbH